

## 西宝生物——储能材料之锂电池



锂电池是由正极材料、负极材料、电解质和隔膜组成。其中，正极材料和负极材料是锂电池的关键材料。

一般锂电池可以分为两大类，分别是锂金属电池和锂离子电池。这是根据锂存在的形态来定义的，锂金属电池是用金属锂做电极，而锂离子电池则是以离子形态存在于电极。

新型锂电池材料主要包括以下几种原材料：

1. 正极材料：主要有三元材料（如 $\text{LiCoO}_2$ 、 $\text{LiNiCoAlO}_2$ 等）、钴酸锂、锰酸锂、铁酸锂、磷酸铁锂等。
2. 负极材料：主要有石墨、硅、锡、碳纳米管等。
3. 电解质：主要有有机电解质（如聚合物电解质、液态电解质等）和无机电解质（如氧化物电解质等）。
4. 隔膜：主要有聚丙烯膜、聚酰亚胺膜等。

锂金属电池是以锂金属或锂合金为负极材料，采用非水电解质溶液的电池。锂金属电池是一次性电池，通过金属锂的腐蚀或氧化来产生电能，用完就废了，不能充电。

锂离子电池是可充电电池，正极通常采用聚偏二氟乙烯做粘结剂，溶剂是N-甲基吡咯烷酮，负极采用非水的有机溶剂做电解液。主要依靠锂离子在正极和负极之间移动来工作。在充放电过程中， $\text{Li}^+$ 在两个电极之间往返嵌入和脱嵌：充电时， $\text{Li}^+$ 从正极脱嵌，经过电解质嵌入负极，负极处于富锂状态；放电时则相反。

需要注意的是，锂电池材料的选择需要根据电池的使用场景、性能要求、成本等因素来综合考虑，不同的材料组合会影响电池的性能和成本。

西宝生物提供储能电池原材料，欢迎来电咨询：400-021-8158 / 021-50272975。

中文名称	英文名称	CAS
3,4-乙烯二氧噻吩(EDOT)	3,4-Ethylenedioxythiophene	126213-50-1
4-苯基咪唑	4-Phenylimidazole	670-95-1
Boc-L-β-高精氨酸对甲苯磺酸盐	BOC-BETA-HOMOARG(TOS)-OH	136271-81-3
对甲苯磺酸铁	Iron(III) p-toluenesulfonate	77214-82-5
硫酸锂一水合物	Lithium sulfate monohydrate	10102-25-7
硫酸锰一水合物	Manganese(II) sulfate monohydrate	10034-96-5
甲基乙烯基醚/马来酸酐共聚物	Poly(methyl vinyl ether-alt-maleic anhydride)	9011-16-9
聚苯乙烯磺酸钠 (PSS)	Poly(sodium-p-styrenesulfonate)	25704-18-1
聚乙烯磺酸-苯乙烯磺酸	PVS-SS	
S-EDOT	sodium 3-(2,3-dihydrothieno[3,4-b][1,4]dioxin-2-yl)methoxy-1-methylpropanesulfonate	1526907-40-3
(E,E,E,E)-1,2,4,5-四[4-[双(4-甲氧基苯基)氨基]苯乙炔基]苯	(E,E,E,E)-1,2,4,5-Tetrakis[4-[bis(4-methoxyphenyl)amino]styryl]benzene	
(二茂铁甲基)三甲基氯化铵	(Ferrocenylmethyl)trimethylammonium Chloride	83617-79-2
[4-(3,6-二甲氧基-9H-咔唑-9-基)丁基]膦酸	[4-(3,6-Dimethoxy-9H-carbazol-9-yl)butyl]phosphonic Acid	2922526-56-3
[4-(3,6-二甲基-9H-咔唑-9-基)丁基]膦酸	[4-(3,6-Dimethyl-9H-carbazol-9-yl)butyl]phosphonic Acid	2747959-96-0
1-萘甲基溴化铵	1-Naphthylmethylammonium Bromide	217309-83-6
1-(2-羟基乙基)-3-甲基咪唑双(三氟甲磺酰基)亚胺	1-(2-Hydroxyethyl)-3-methylimidazolium Bis(trifluoromethanesulfonyl)imide	174899-86-6
1-(2-羟乙基)-3-甲基氯化咪唑	1-(2-Hydroxyethyl)-3-methylimidazolium Chloride	61755-34-8
1-(2-羟乙基)-3-甲基咪唑四氟硼酸盐	1-(2-Hydroxyethyl)-3-methylimidazolium Tetrafluoroborate	374564-83-7
1-(2-甲氧乙基)-1-甲基吡咯烷鎓双(氟磺酰基)亚胺	1-(2-Methoxyethyl)-1-methylpyrrolidinium Bis(fluorosulfonyl)imide	1235234-47-5
1,1,1-三氟-N,N-二甲基甲磺酰胺	1,1,1-Trifluoro-N,N-dimethylmethanesulfonamide	28048-17-1
1,12-十二胺氢溴酸盐	1,12-Dodecanediammonium Dibromide	185671-15-2
1,12-十二胺氢碘酸盐	1,12-Dodecanediammonium Diiodide	342888-94-2
1,1'-双[3-(三甲基铵)丙基]-4,4'-联吡啶四氯化物二水合物	1,1'-Bis[3-(trimethylammonio)propyl]-4,4'-bipyridinium Tetrachloride Dihydrate	108228-37-1
1,1'-双[3-(三甲基铵)丙基]二茂铁二氯化物	1,1'-Bis[3-(trimethylammonio)propyl]ferrocene Dichloride	2093414-16-3
1,1-二苯基-2-苦基肼自由基	1,1-Diphenyl-2-picrylhydrazyl Free Radical	1898-66-4
1,2,2,3-丙烷四甲腈	1,2,2,3-Propanetetracarbonitrile	1274904-48-1
1,2-环己二胺氢溴酸盐	1,2-Cyclohexanediammonium dibromide	
1,2-环己二胺氢碘酸盐	1,2-cyclohexanediammonium diiodide	
1,2-二甲基-3-丙基碘化咪唑鎓	1,2-Dimethyl-3-propylimidazolium Iodide	218151-78-1
1,2-苯二胺氢碘酸盐	1,2-Phenyldiammonium diiodide	81584-89-6
1,3-环己二胺氢碘酸盐	1,3-Cyclohexyldimethylammonium diiodide	
1,3,2-二氧硫杂环戊烷2,2-二氧化物	1,3,2-Dioxathiolane 2,2-Dioxide	1072-53-3
1,3,2-二氧硫杂戊烷2-氧化物	1,3,2-Dioxathiolane 2-Oxide	3741-38-6
1,3,5-三[4-[双(4-甲氧基苯基)氨基]苯基]苯	1,3,5-Tris[4-[bis(4-methoxyphenyl)amino]phenyl]benzene	142894-38-0
1,3,6-己烷三腈	1,3,6-Hexanetricarbonitrile	1772-25-4
1,3-环己二胺氢溴酸盐	1,3-Cyclohexyldimethylammonium dibromide	
1,3-二氨基丙烷二氢溴酸盐	1,3-Diaminopropane Dihydrobromide	18773-03-0
1,3-二氨基丙烷二盐酸盐	1,3-Diaminopropane Dihydrochloride	10517-44-9
1,3-二氨基丙烷二碘酸盐	1,3-Diaminopropane Dihydroiodide	120675-53-8
1,3-二甲基-1H-咪唑-3-鎓四氟硼酸盐	1,3-Dimethyl-1H-imidazol-3-ium Tetrafluoroborate	121091-31-4

中文名称	英文名称	CAS
1,3-二甲基咪唑双(三氟甲磺酰)亚胺	1,3-Dimethylimidazolium Bis(trifluoromethanesulfonyl)imide	174899-81-1
氯化1,3-二甲基咪唑	1,3-Dimethylimidazolium Chloride	79917-88-7
1,3-二甲基咪唑鎓二甲磷酸酯	1,3-Dimethylimidazolium Dimethyl Phosphate	654058-04-5
1,3-二甲基咪唑鎓碘化物	1,3-Dimethylimidazolium Iodide	4333-62-4
1,3-二甲基咪唑鎓甲基硫酸盐	1,3-Dimethylimidazolium Methyl Sulfate	97345-90-9
1,3-苯二胺氢碘酸盐	1,3-Phenyldiammonium diiodide	
1,3-苯二甲胺氢碘酸盐	1,3-Phenyldimethylammonium diiodide	
1,3-丙烷磺内酯	1,3-Propanesultone	1120-71-4
1,4-环己二胺氢碘酸盐	1,4-Cyclohexyldimethylammonium diiodide	
1,4-丁磺酸内酯	1,4-Butanesultone	1633-83-6
1,4-环己二胺氢溴酸盐	1,4-Cyclohexyldimethylammonium dibromide	
1,4-丁二胺氢碘酸盐	1,4-Diaminobutane Dihydroiodide	916849-52-0
1,4-二氮杂双环[2.2.2]辛烷二氢碘酸盐	1,4-Diazabicyclo[2.2.2]octane Dihydriodide	33322-06-4
1,4-二氮杂双环[2.2.2]辛烷二氢溴酸盐	1,4-Diazabicyclo[2.2.2]octane Dihydrobromide	54581-69-0
1,4-二氮杂双环[2.2.2]辛烷二盐酸盐	1,4-Diazabicyclo[2.2.2]octane Dihydrochloride	49563-87-3
1,4-二-叔-丁基-2,5-双(2-甲氧基乙氧基)苯	1,4-Di-tert-butyl-2,5-bis(2-methoxyethoxy)benzene	1350770-63-6
1,4-苯二胺二氢碘化物	1,4-Phenylenediamine Dihydriodide	116469-02-4
1,8-辛二胺盐酸盐	1,8-Octanediammonium dichloride	7613-16-3
10-[2-(2-甲氧基乙氧基)乙基]-10H-吩噻嗪	10-[2-(2-Methoxyethoxy)ethyl]-10H-phenothiazine	2098786-35-5
10-乙酰苯噻嗪	10-Acetylphenothiazine	1628-29-1
1-烯丙基-1-甲基吡咯烷鎓双(三氟甲磺酰)亚胺	1-Allyl-1-methylpyrrolidinium Bis(trifluoromethanesulfonyl)imide	1059624-23-5
1-烯丙基-3-甲基-1H-咪唑-3-鎓溴化物	1-Allyl-3-methyl-1H-imidazol-3-ium Bromide	31410-07-8
1-烯丙基-3-甲基咪唑双(三氟甲磺酰)亚胺	1-Allyl-3-methylimidazolium Bis(trifluoromethanesulfonyl)imide	655249-87-9
1-烯丙基-3-甲基氯化咪唑	1-Allyl-3-methylimidazolium Chloride	65039-10-3
1-苄基-3-甲基咪唑氯化物	1-Benzyl-3-methylimidazolium Chloride	36443-80-8
1-苄基-3-甲基咪唑六氟磷酸盐	1-Benzyl-3-methylimidazolium	433337-11-2
1-苄基-3-甲基咪唑四氟硼酸盐	1-Benzyl-3-methylimidazolium Tetrafluoroborate	500996-04-3
1-丁基-1-甲基哌啶双(三氟甲磺酰)亚胺盐	1-Butyl-1-methylpiperidinium Bis(trifluoromethanesulfonyl)imide	623580-02-9
1-丁基-1-甲基溴化哌啶鎓	1-Butyl-1-methylpiperidinium Bromide	94280-72-5
1-丁基-1-甲基吡咯烷鎓双(氟磺酰)亚胺	1-Butyl-1-methylpyrrolidinium	1057745-51-3
1-丁基-1-甲基吡咯烷鎓双(三氟甲基磺酰基)亚胺	1-Butyl-1-methylpyrrolidinium bis(trifluoromethylsulfonyl)imide	223437-11-4
1-丁基-1-甲基吡咯烷鎓溴化物	1-Butyl-1-methylpyrrolidinium Bromide	93457-69-3
1-丁基-1-甲基吡咯烷氯化物	1-Butyl-1-methylpyrrolidinium Chloride	479500-35-1
1-丁基-1-甲基吡咯烷二氰胺盐	1-Butyl-1-methylpyrrolidinium Dicyanamide	370865-80-8
1-丁基-1-甲基吡咯烷六氟磷酸盐	1-Butyl-1-methylpyrrolidinium	330671-29-9
1-丁基-1-甲基-1-吡咯烷鎓四氟硼酸盐	1-Butyl-1-methylpyrrolidinium tetrafluoroborate	345984-11-4
1-丁基-1-甲基吡咯三氟甲磺酸盐	1-Butyl-1-methylpyrrolidinium Trifluoromethanesulfonate	367522-96-1
1-丁基-2,3-二甲基咪唑双(三氟甲基磺酰)亚胺	1-Butyl-2,3-dimethylimidazolium Bis(trifluoromethanesulfonyl)imide	350493-08-2
1-丁基-2,3-二甲基氯化咪唑鎓	1-Butyl-2,3-dimethylimidazolium Chloride	98892-75-2
1-丁基-2,3-二甲基咪唑鎓六氟磷酸盐	1-Butyl-2,3-dimethylimidazolium	227617-70-1
1-丁基-2,3-二甲基咪唑鎓四氟硼酸盐	1-Butyl-2,3-dimethylimidazolium	402846-78-0
1-丁基-2,3-二甲基咪唑鎓三氟甲烷磺酸盐	1-Butyl-2,3-dimethylimidazolium Trifluoromethanesulfonate	765910-73-4

中文名称	英文名称	CAS
1-丁基-3-甲基咪唑乙酸盐	1-Butyl-3-methylimidazolium Acetate	284049-75-8
1-丁基-3-甲基咪唑双(三氟甲磺酰)亚胺	1-Butyl-3-methylimidazolium Bis(trifluoromethanesulfonyl)imide	174899-83-3
1-丁基-3-甲基溴化咪唑鎓	1-Butyl-3-methylimidazolium Bromide	85100-77-2
1-丁基-3-甲基氯化咪唑鎓 (BMIMCl)	1-Butyl-3-methylimidazolium chloride	79917-90-1
1-丁基-3-甲基咪唑二氰胺盐	1-Butyl-3-methylimidazolium Dicyanamide	448245-52-1
1-丁基-3-甲基咪唑六氟锑酸盐	1-Butyl-3-methylimidazolium	174645-81-9
1-丁基-3-甲基六氟磷酸盐咪唑鎓	1-Butyl-3-methylimidazolium	174501-64-5
1-丁基-3-甲基咪唑硫酸氢盐	1-Butyl-3-methylimidazolium Hydrogen Sulfate	262297-13-2
1-丁基-3-甲基碘化咪唑鎓	1-Butyl-3-methylimidazolium Iodide	65039-05-6
1-丁基-3-甲基咪唑甲磺酸盐	1-Butyl-3-methylimidazolium Methanesulfonate	342789-81-5
1-丁基-3-甲基咪唑硫酸甲酯	1-Butyl-3-methylimidazolium Methyl Sulfate	401788-98-5
1-丁基-3-甲基咪唑鎓四氯高铁酸盐	1-Butyl-3-methylimidazolium Tetrachloroferrate	359845-21-9
1-丁基-3-甲基咪唑四氟硼酸盐	1-Butyl-3-methylimidazolium Tetrafluoroborate	174501-65-6
1-丁基-3-甲基咪唑鎓硫氰酸盐	1-Butyl-3-methylimidazolium Thiocyanate	344790-87-0
1-丁基-3-甲基咪唑鎓三溴化物	1-Butyl-3-methylimidazolium Tribromide	820965-08-0
1-丁基-3-甲基咪唑鎓三氰基甲烷化物	1-Butyl-3-methylimidazolium Tricyanomethanide	878027-73-7
1-丁基-3-甲基咪唑啉三氟(三氟甲基)硼酸盐	1-Butyl-3-methylimidazolium Trifluoro(trifluoromethyl)borate	741677-68-9
1-丁基-3-甲基咪唑三氟乙酸盐	1-Butyl-3-methylimidazolium Trifluoroacetate	174899-94-6
1-丁基-3-甲基咪唑三氟甲烷磺酸盐	1-Butyl-3-methylimidazolium	174899-66-2
1-丁基-3-甲基吡啶溴化物	1-Butyl-3-methylpyridinium Bromide	26576-85-2
1-丁基-3-甲基氯化吡啶鎓	1-Butyl-3-methylpyridinium Chloride	125652-55-3
1-丁基-3-乙烯基-1H-咪唑-3-鎓溴化物	1-Butyl-3-vinyl-1H-imidazol-3-ium Bromide	34311-90-5
1-丁基-4-甲基吡啶双(三氟甲磺酰)亚胺	1-Butyl-4-methylpyridinium Bis(trifluoromethanesulfonyl)imide	475681-62-0
1-丁基-4-甲基吡啶溴盐	1-Butyl-4-methylpyridinium Bromide	65350-59-6
1-丁基-4-甲基氯化吡啶	1-Butyl-4-methylpyridinium Chloride	112400-86-9
1-丁基-4-甲基吡啶六氟磷酸盐	1-Butyl-4-methylpyridinium Hexafluorophosphate	401788-99-6
1-丁基-4-甲基吡啶四氟硼酸盐	1-Butyl-4-methylpyridinium Tetrafluoroborate	343952-33-0
1-丁基吡啶双(三氟甲烷磺酰)亚胺盐	1-Butylpyridinium	187863-42-9
1-丁基溴化吡啶	1-Butylpyridinium Bromide	874-80-6
1-丁基氯化吡啶	1-Butylpyridinium Chloride	1124-64-7
1-丁基吡啶六氟磷酸盐	1-Butylpyridinium Hexafluorophosphate	186088-50-6
1-丁基吡啶四氟硼酸盐	1-Butylpyridinium Tetrafluoroborate	203389-28-0
1-癸基-3-甲基咪唑双(三氟甲烷磺酰基)酰亚胺	1-Decyl-3-methylimidazolium Bis(trifluoromethanesulfonyl)imide	433337-23-6
1-癸基-3-甲基溴化咪唑	1-Decyl-3-methylimidazolium Bromide	188589-32-4
1-癸基-3-甲基氯化咪唑	1-Decyl-3-methylimidazolium Chloride	171058-18-7
1-癸基-3-甲基咪唑四氟硼酸盐	1-Decyl-3-methylimidazolium Tetrafluoroborate	244193-56-4
1-十二烷基-3-甲基-1H-咪唑-3-鎓氯化物	1-Dodecyl-3-methyl-1H-imidazol-3-ium Chloride	114569-84-5
1-十二烷基-3-甲基咪唑鎓双(三氟甲磺酰基)亚胺	1-Dodecyl-3-methylimidazolium Bis(trifluoromethanesulfonyl)imide	404001-48-5
1-十二烷基-3-甲基咪唑溴化物	1-Dodecyl-3-methylimidazolium Bromide	61546-00-7
1-乙基-1-甲基吡咯烷鎓双(三氟甲基磺酰基)亚胺	1-Ethyl-1-methylpyrrolidinium bis(trifluoromethylsulfonyl)imide	223436-99-5
1-乙基-1-甲基溴化吡咯烷	1-Ethyl-1-methylpyrrolidinium Bromide	69227-51-6
1-乙基-1-甲基吡咯烷鎓四氟硼酸盐	1-Ethyl-1-methylpyrrolidinium Tetrafluoroborate	117947-85-0
1-乙基-2,3-二甲基咪唑鎓双(三氟甲磺酰)亚胺	1-Ethyl-2,3-dimethylimidazolium Bis(trifluoromethanesulfonyl)imide	174899-90-2
1-乙基-2-甲基溴化吡啶	1-Ethyl-2-methylpyridinium Bromide	32353-50-7

中文名称	英文名称	CAS
1-乙基-3-(羟甲基)吡啶鎓乙磺酸盐	1-Ethyl-3-(hydroxymethyl)pyridinium Ethyl Sulfate	1064704-03-5
1-乙基-3-甲基咪唑鎓乙酸盐	1-Ethyl-3-methylimidazolium Acetate	143314-17-4
1-乙基-3-甲基咪唑双三氟甲磺酰亚胺盐	1-Ethyl-3-methylimidazolium bis(trifluoromethylsulfonyl)imide	174899-82-2
1-乙基-3-甲基溴化咪唑鎓	1-Ethyl-3-methylimidazolium Bromide	65039-08-9
1-乙基-3-甲基氯化咪唑鎓	1-Ethyl-3-methylimidazolium Chloride	65039-09-0
1-乙基-3-甲基咪唑鎓二氰胺盐	1-Ethyl-3-methylimidazolium Dicyanamide	370865-89-7
1-乙基-3-甲基咪唑二乙基磷酸盐	1-Ethyl-3-methylimidazolium Diethyl Phosphate	848641-69-0
1-乙基-3-甲基咪唑鎓乙基磺酸盐	1-Ethyl-3-methylimidazolium Ethyl Sulfate	342573-75-5
1-乙基-3-甲基咪唑六氟磷酸盐	1-Ethyl-3-methylimidazolium Hexafluorophosphate	155371-19-0
1-乙基-3-甲基碘化咪唑鎓	1-Ethyl-3-methylimidazolium Iodide	35935-34-3
1-乙基-3-甲基咪唑鎓甲烷磺酸盐	1-Ethyl-3-methylimidazolium Methanesulfonate	145022-45-3
1-乙基-3-甲基咪唑甲基硫酸盐	1-Ethyl-3-methylimidazolium Methyl Sulfate	516474-01-4
1-乙基-3-甲基咪唑硝酸盐	1-Ethyl-3-methylimidazolium Nitrate	143314-14-1
1-乙基-3-甲基咪唑鎓对甲苯磺酸盐	1-Ethyl-3-methylimidazolium p-Toluenesulfonate	328090-25-1
1-乙基-3-甲基咪唑鎓四氯铝酸盐	1-Ethyl-3-methylimidazolium tetrachloroaluminate	80432-05-9
1-乙基-3-甲基咪唑鎓四氯高铁酸盐	1-Ethyl-3-methylimidazolium Tetrachloroferrate	850331-04-3
1-乙基-3-甲基咪唑四氟硼酸盐	1-Ethyl-3-methylimidazolium tetrafluoroborate	143314-16-3
1-乙基-3-甲基咪唑硫氰酸盐	1-Ethyl-3-methylimidazolium Thiocyanate	331717-63-6
1-乙基-3-甲基咪唑鎓三氰基甲烷化物	1-Ethyl-3-methylimidazolium Tricyanomethanide	666823-18-3
1-乙基-3-甲基咪唑啉三氟(三氟甲基)硼酸盐	1-Ethyl-3-methylimidazolium Trifluoro(trifluoromethyl)borate	681856-28-0
1-乙基-3-甲基咪唑三氟乙酸盐	1-Ethyl-3-methylimidazolium Trifluoroacetate	174899-65-1
1-乙基-3-甲基咪唑三氟甲烷磺酸盐	1-Ethyl-3-methylimidazolium	145022-44-2
1-乙基-3-甲基吡啶鎓双(三氟甲磺酰)亚胺	1-Ethyl-3-methylpyridinium Bis(trifluoromethanesulfonyl)imide	841251-37-4
1-乙基-3-甲基吡啶鎓乙磺酸盐	1-Ethyl-3-methylpyridinium Ethyl Sulfate	872672-50-9
1-乙基-4-甲基溴化吡啶	1-Ethyl-4-methylpyridinium Bromide	32353-49-4
1-乙基溴化吡啶鎓	1-Ethylpyridinium Bromide	1906-79-2
1-乙基氯化吡啶鎓	1-Ethylpyridinium Chloride	2294-38-4
庚基溴化铵	1-Heptan ammonium Bromide	15590-85-9
庚基碘化铵	1-Heptan ammonium Iodide	60734-62-5
1-十六烷基-3-甲基氯化咪唑	1-Hexadecyl-3-methylimidazolium Chloride	61546-01-8
1-己基-1,4-二氮杂双环[2.2.2]辛烷-1-鎓碘化物	1-Hexyl-1,4-diazabicyclo[2.2.2]octan-1-ium Iodide	1009321-13-4
1-己基-1-甲基吡咯烷-1-鎓双((三氟甲基)磺酰基)亚胺	1-Hexyl-1-methylpyrrolidin-1-ium Bis((trifluoromethyl)sulfonyl)amide	380497-19-8
1-己基-2,3-二甲基咪唑啉碘化物	1-Hexyl-2,3-dimethylimidazolium Iodide	288627-94-1
1-己基-3-甲基咪唑双(三氟甲烷磺酰基)酰亚胺	1-Hexyl-3-methylimidazolium Bis(trifluoromethanesulfonyl)imide	382150-50-7
1-己基-3-甲基咪唑鎓溴化物	1-Hexyl-3-methylimidazolium Bromide	85100-78-3
1-己基-3-甲基氯化咪唑鎓	1-Hexyl-3-methylimidazolium Chloride	171058-17-6
1-己基-3-甲基咪唑鎓六氟磷酸盐	1-Hexyl-3-methylimidazolium	304680-35-1
1-己基-3-甲基咪唑碘化物	1-Hexyl-3-methylimidazolium Iodide	178631-05-5
1-己基-3-甲基咪唑鎓四氟硼酸盐	1-Hexyl-3-methylimidazolium Tetrafluoroborate	244193-50-8
1-己基-3-甲基咪唑鎓三氟甲烷磺酸盐	1-Hexyl-3-methylimidazolium Trifluoromethanesulfonate	460345-16-8
1-己基吡啶鎓六氟磷酸盐	1-Hexylpyridinium Hexafluorophosphate	797789-00-5
1-甲基-1H-咪唑-3-鎓六氟磷酸盐	1-Methyl-1H-imidazol-3-ium Hexafluorophosphate	57367-08-5
1-甲基-1-n-辛基吡咯烷鎓双(三氟甲磺酰基)亚胺	1-Methyl-1-n-octylpyrrolidinium Bis(trifluoromethanesulfonyl)imide	927021-43-0
1-甲基-1-辛基吡咯烷-1-鎓氯化物	1-Methyl-1-octylpyrrolidin-1-ium Chloride	909398-60-3

中文名称	英文名称	CAS
1-甲基-1-戊基吡咯烷鎓双(三氟甲磺酰基)亚胺	1-Methyl-1-pentylpyrrolidinium Bis(trifluoromethanesulfonyl)imide	380497-17-6
1-甲基-1-丙基哌啶鎓双(氟磺酰)亚胺	1-Methyl-1-propylpiperidinium	911303-46-3
1-甲基-1-丙基哌啶-1-鎓双((三氟甲基)磺酰基)酰亚胺	1-Methyl-1-propylpiperidinium bis(trifluoromethylsulfonyl)imide	608140-12-1
1-甲基-1-丙基哌啶鎓溴化物	1-Methyl-1-propylpiperidinium Bromide	88840-42-0
1-甲基-1-丙基吡咯烷-1-鎓溴化物	1-Methyl-1-propylpyrrolidin-1-ium Bromide	608140-09-6
1-甲基-1-丙基吡咯烷鎓双(氟磺酰)亚胺	1-Methyl-1-propylpyrrolidinium Bis(fluorosulfonyl)imide	852620-97-4
1-甲基-1-丙基吡咯烷双(三氟甲磺酰基)亚胺	1-Methyl-1-propylpyrrolidinium Bis(trifluoromethanesulfonyl)imide	223437-05-6
1-甲基-3-(4-磺丁基)咪唑鎓双(三氟甲磺酰基)亚胺	1-Methyl-3-(4-sulfobutyl)imidazolium Bis(trifluoromethanesulfonyl)imide	909390-59-6
1-甲基-3-(4-磺丁基)咪唑硫酸氢盐	1-Methyl-3-(4-sulfobutyl)imidazolium Hydrogen	827320-59-2
1-甲基-3-[3-(三甲氧基硅基)丙基]-1H-咪唑-3-鎓氯化物	1-Methyl-3-[3-(trimethoxysilyl)propyl]-1H-imidazol-3-ium Chloride	856925-70-7
1-甲基-3-正辛基咪唑双(三氟甲磺酰基)亚胺	1-Methyl-3-n-octylimidazolium Bis(trifluoromethanesulfonyl)imide	178631-04-4
1-甲基-3-正辛基咪唑鎓溴化物	1-Methyl-3-n-octylimidazolium Bromide	61545-99-1
1-甲基-3-正辛基咪唑鎓氯化物	1-Methyl-3-n-octylimidazolium Chloride	64697-40-1
1-甲基-3-正辛基咪唑鎓六氟磷酸盐	1-Methyl-3-n-octylimidazolium	304680-36-2
1-甲基-3-正辛基咪唑四氟硼酸盐	1-Methyl-3-n-octylimidazolium Tetrafluoroborate	244193-52-0
1-甲基-3-正辛基咪唑鎓三氟甲烷磺酸盐	1-Methyl-3-n-octylimidazolium Trifluoromethanesulfonate	403842-84-2
1-甲基-3-戊基溴化咪唑	1-Methyl-3-pentylimidazolium Bromide	343851-31-0
1-甲基-3-苯基丙胺溴	1-Methyl-3-phenylpropylammonium Bromide	
1-甲基-3-苯基丙胺碘化铵	1-Methyl-3-phenylpropylammonium Iodide	
1-甲基-3-丙基咪唑双(三氟甲磺酰基)亚胺	1-Methyl-3-propylimidazolium Bis(trifluoromethanesulfonyl)imide	216299-72-8
1-甲基-3-丙基溴化咪唑	1-Methyl-3-propylimidazolium Bromide	85100-76-1
1-甲基-3-丙基氯化咪唑 (PMIMCl)	1-Methyl-3-propylimidazolium Chloride	79917-89-8
1-甲基-3-丙基碘化咪唑鎓	1-Methyl-3-propylimidazolium Iodide	119171-18-5
1-甲基-3-丙基咪唑四氟硼酸盐	1-Methyl-3-propylimidazolium Tetrafluoroborate	244193-48-4
1-甲基咪唑双(三氟甲磺酰基)亚胺	1-Methylimidazole	353239-08-4
1-甲基咪唑氢溴酸盐	1-Methylimidazole Hydrobromide	101023-58-9
1-甲基咪唑三氟甲磺酸盐	1-Methylimidazole Trifluoromethanesulfonate	99257-94-0
1-甲基吡啶鎓双(三氟甲磺酰基)亚胺	1-Methylpyridinium	742086-10-8
1-甲基吡啶鎓六氟磷酸盐	1-Methylpyridinium Hexafluorophosphate	117029-34-2
1-萘甲基碘化铵	1-Naphthylmethylammonium iodide	879678-26-9
1-戊胺氢溴酸盐	1-Pentanamine Hydrobromide	7334-94-3
1-戊胺盐酸盐	1-Pentanamine Hydrochloride	142-65-4
1-丙基氯化吡啶	1-Propylpyridinium Chloride	23271-47-8
1-乙烯基咪唑双(三氟甲磺酰基)亚胺	1-Vinylimidazole	1013027-27-4
2-萘甲基溴化铵	2-Naphthylmethylammonium Bromide	
2-萘甲基碘化铵	2-Naphthylmethylammonium iodide	
2-(1-吡唑基)吡啶	2-(1-Pyrazolyl)pyridine	25700-11-2
2-(1-吡咯烷基)乙胺二盐酸盐	2-(1-Pyrrolidinyl)ethanamine Dihydrochloride	65592-36-1
2-(2,2,2-三氟乙氧基)-1,3,2-二氧磷杂环戊烷 2-氯化物	2-(2,2,2-Trifluoroethoxy)-1,3,2-dioxaphospholane 2-Oxide	67605-68-9
2-(4-氟苯基)乙胺氢溴酸盐	2-(4-Fluorophenyl)ethylamine Hydrobromide	1807536-06-6
2-(4-氟苯基)乙胺碘酸盐	2-(4-Fluorophenyl)ethylamine Hydroiodide	1413269-55-2

中文名称	英文名称	CAS
2-(4-甲氧基苯基)乙胺氢溴酸盐	2-(4-Methoxyphenyl)ethylamine Hydrobromide	55-81-2 (free base)
2-(4-甲氧基苯基)乙胺碘酸盐	2-(4-Methoxyphenyl)ethylamine Hydroiodide	
2,2'-(1,3-苯基)二[5-(4-叔丁基苯基)-1,3,4-噁二唑]	2,2'-(1,3-Phenylene)bis[5-(4-tert-butylphenyl)-1,3,4-oxadiazole]	138372-67-5
2,2'-(2-乙烯基蒽-9,10-二亚基)双(1,3-二硫杂环戊二烯)	2,2'-(2-Vinylanthracene-9,10-diylidene)bis(1,3-dithiole)	736998-56-4
2,2'-(2-乙烯基蒽-9,10-二亚基)二丙二腈	2,2'-(2-Vinylanthracene-9,10-	1612793-07-3
2,2,2-三氟乙胺氢溴酸盐	2,2,2-Trifluoroethylammonium Bromide	373-87-5
2,2,2-三氟乙胺碘酸盐	2,2,2-Trifluoroethylammonium Iodide	1650572-03-4
2,2,6,6-四甲基-4-(2-丙炔基氧)哌啶1-氧自由基	2,2,6,6-Tetramethyl-4-(2-propynoxy)piperidine 1-Oxyl Free Radical	147045-24-7
2,2,6,6-四甲基哌啶1-氧自由基 (升华提纯)	2,2,6,6-Tetramethylpiperidine 1-Oxyl Free Radical (purified by sublimation)	2564-83-2
2,2',7,7'-四[N-4-甲氧基苯基氨基]-2,2',7,7'-四[N-(3-氟-4-甲氧基苯基)]-9,9'-螺二芴	2,2',7,7'-Tetrakis[N,N-4-methoxyphenylamino]-2,2',7,7'-Tetrakis[N,N-3-fluoro-4-methoxyphenyl]-9,9'-spirobifluorene	2644017-35-4
2,2'-(乙烯二氧)双乙胺氢溴酸盐	2,2'-(Ethylenedioxy)diethylammonium dibromide	
2,2'-(乙烯二氧)双乙胺盐酸盐	2,2'-(Ethylenedioxy)diethylammonium dichloride	80685-25-2
2,2'-二喹啉-4,4'-二羧酸	2,2'-Biquinoline-4,4'-dicarboxylic acid	1245-13-2
2,3,4,5,6-五(9-咪唑基)-苯腈	2,3,4,5,6-tetrakis(Carbazol-9-yl)-cyanobenzene	1469700-24-0
2,3-二甲基-1-丙基咪唑鎓双(三氟甲磺酰基)酰亚胺	2,3-Dimethyl-1-propylimidazolium Bis(trifluoromethanesulfonyl)imide	169051-76-7
2,4,4-三甲基戊-2-胺氢溴酸盐	2,4,4-Trimethylpentan-2-amine Hydrobromide	1093859-61-0
2,4,4-三甲基戊-2-胺盐酸盐	2,4,4-Trimethylpentan-2-amine Hydrochloride	58618-91-0
2,4,4-三甲基戊-2-胺碘酸盐	2,4,4-Trimethylpentan-2-amine Hydroiodide	
2,4,6-三甲氧基硼氧六环	2,4,6-Trimethoxyboroxin	102-24-9
2,5-双(4,4'-二甲氧基三苯胺基)-3,4-乙烯二氧噻吩	2,5-bis(4,4'-bis(methoxyphenyl)aminophen-4-yl)-3,4-ethylenedioxythiophene	1622008-73-4
2,6-DBEAQ	2,6-DBEAQ	2370885-20-2
2,6-DPPEAQ	2,6-DPPEAQ	2370885-23-5
噻吩-2-甲脒碘酸盐	2-Amidinothiophene hydroiodide	2568582-25-0
磷酸2-乙基己基二苯酯	2-Ethylhexyl Diphenyl Phosphate	1241-94-7
2-氟乙胺盐酸盐	2-Fluoroethylamine hydrochloride	460-08-2
2-氟苯乙基溴化胺	2-Fluorophenylethylammonium Bromide	
2-氟苯乙基碘化胺	2-Fluorophenylethylammonium Iodide	
2-甲氧基乙胺氢溴酸盐	2-Methoxyethylamine Hydrobromide	663941-77-3
2-甲氧基苯乙基溴化胺	2-Methoxyphenylethylammonium Bromide	
2-甲氧基苯乙基碘化胺	2-Methoxyphenylethylammonium Iodide	
[2-(9H-咪唑-9-基)乙基]膦酸	2PACz	20999-38-6
2-苯乙胺碘酸盐	2-Phenylethylamine Hydroiodide	151059-43-7
2-炔丙胺碘酸盐	2-Propynylammonium iodide	205817-86-3
2-噻吩甲胺氢溴酸盐	2-Thiophenemethylammonium Bromide	59852-33-4
2-噻吩甲胺碘酸盐	2-Thiophenemethylammonium Iodide	2247100-32-7
2-噻吩乙胺氢溴酸盐	2-Thiophenethylammonium Bromide	2490324-60-0
2-噻吩乙胺碘酸盐	2-Thiophenethylammonium Iodide	2414055-94-8
2-乙烯基蒽醌	2-Vinylanthraquinone	13388-33-5
3-(羧甲基)-1-甲基-1H-咪唑-3-鎓氯化物	3-(Carboxymethyl)-1-methyl-1H-imidazol-3-ium Chloride	700370-07-6
3,3'-(丁烷-1,4-二基)双(1-乙烯基-3-咪唑鎓)双(三氟甲磺酰)亚胺	3,3'-(Butane-1,4-diyl)bis(1-vinyl-3-imidazolium) Bis(trifluoromethanesulfonyl)imide	1312310-16-9

中文名称	英文名称	CAS
3,3'-氧二丙腈	3,3'-Oxydipropionitrile	1656-48-0
3,3-二苯基丙胺氢碘酸盐	3,3-Diphenylpropylammonium Iodide	2494199-73-2
3-氨基己酸氢溴酸盐	3-Aminohexanoic Acid Hydrobromide	
3-氨基己酸氢碘酸盐	3-Aminohexanoic Acid Hydroiodide	
3-氨基丙酸氢碘酸盐	3-Aminoprotanoic Acid Hydroiodide	
3-丁基-1-甲基-1H-咪唑-3-鎓甲基苯磺酸盐	3-Butyl-1-methyl-1H-imidazol-3-ium 4-Methylbenzenesulfonate	410522-18-8
3-羧基-2,2,5,5-四甲基吡咯烷-1-氧基 自由基	3-Carboxy-2,2,5,5-tetramethylpyrrolidine 1-Oxyl Free Radical	2154-68-9
3-乙氧基丙腈	3-Ethoxypropionitrile	2141-62-0
3-乙基-1-甲基-1H-咪唑-3-鎓双(氟磺酰基)亚胺盐	3-Ethyl-1-methyl-1H-imidazol-3-ium Bis(fluorosulfonyl)azanide	235789-75-0
3-乙基-1-甲基-1H-咪唑-3-鎓二甲基磷酸盐	3-Ethyl-1-methyl-1H-imidazol-3-ium Dimethyl	945611-27-8
3-乙基-1-乙烯基-1H-咪唑-3-鎓溴化物	3-Ethyl-1-vinyl-1H-imidazol-3-ium Bromide	34311-88-1
3-乙基-1-乙烯基咪唑鎓双(三氟甲磺酰)亚胺	3-Ethyl-1-vinylimidazolium Bis(trifluoromethanesulfonyl)imide	319476-28-3
3-氟苯乙基溴化胺	3-Fluorophenylethylammonium Bromide	
3-氟苯乙基氯化胺	3-Fluorophenylethylammonium Chloride	199296-53-2
3-氟苯乙基碘化胺	3-Fluorophenylethylammonium Iodide	
3-十六烷基-1-甲基-1H-咪唑-3-鎓溴化物	3-Hexadecyl-1-methyl-1H-imidazol-3-ium	132361-22-9
3-羟基丙腈	3-Hydroxypropionitrile	109-78-4
3-甲氧基苯乙基溴化胺	3-Methoxyphenylethylammonium Bromide	60189-29-9
3-甲氧基苯乙基碘化胺	3-Methoxyphenylethylammonium Iodide	2760292-10-0
3-甲氧基丙腈	3-Methoxypropionitrile	110-67-8
3-甲基环丁砜	3-Methylsulfolane	872-93-5
3-吡啶甲胺氢溴酸盐	3-PyABr	2144755-94-0
4-(2-氯乙酰氨基)-2,2,6,6-四甲基哌啶-1-氧基 自由基	4-(2-Chloroacetamido)-2,2,6,6-tetramethylpiperidine 1-Oxyl Free Radical	36775-23-2
4-(2-碘代乙酰氨基)-2,2,6,6-四甲基哌啶-1-氧基 自由基	4-(2-Iodoacetamido)-2,2,6,6-tetramethylpiperidine 1-Oxyl Free Radical	25713-24-0
4-(羟甲基)-5-甲基-1,3-二氧戊环-2-酮	4-(Hydroxymethyl)-5-methyl-1,3-dioxol-2-one	91526-18-0
4-(三氟甲基)-1,3-二氧戊环-2-酮	4-(Trifluoromethyl)-1,3-dioxolan-2-one	167951-80-6
4-(三氟甲基)苯胺盐酸盐	4-(Trifluoromethyl)aniline Hydrochloride	90774-69-9
4-(三氟甲基)苄胺氢溴酸盐	4-(Trifluoromethyl)benzylamine Hydrobromide	3300-51-4 (free base)
4-(三氟甲基)苄胺氢碘酸盐	4-(Trifluoromethyl)benzylamine Hydroiodide	
4-乙酰氨基-2,2,6,6-四甲基哌啶-1-氧基 自由基 [氧化反应的催化剂]	4-Acetamido-2,2,6,6-tetramethylpiperidine 1-Oxyl Free Radical [Catalyst for Oxidation]	14691-89-5
4-氨基-2,2,6,6-四甲基哌啶-1-氧基 自由基 (升华提纯)	4-Amino-2,2,6,6-tetramethylpiperidine 1-Oxyl Free Radical (purified by sublimation)	14691-88-4
4-氨基偶氮苯氢溴酸盐	4-Aminoazobenzene hydrobromide	
4-氨基偶氮苯氢碘酸盐	4-Aminoazobenzene hydroiodide	
p-Br-甲基苄胺溴	4-Br-1-Phenylethylammonium Bromide	
p-Br-甲基苄胺碘	4-Br-1-Phenylethylammonium Iodide	
4-溴苯乙基溴化胺	4-Bromophenylethylammonium Bromide	206559-45-7
4-溴苯乙基氯化胺	4-Bromophenylethylammonium Chloride	39260-89-4
4-溴苯乙基碘化胺	4-Bromophenylethylammonium Iodide	2421114-10-3
4-羧基-2,2,6,6-四甲基哌啶-1-氧基 自由基	4-Carboxy-2,2,6,6-tetramethylpiperidine 1-Oxyl Free Radical	37149-18-1
4-氯-1,3-二氧五环-2-酮	4-Chloro-1,3-dioxolan-2-one	3967-54-2
4-氯苯乙基氯化胺	4-Chlorophenylethylammonium Chloride	2492-83-2



中文名称	英文名称	CAS
4-氯苯乙基碘化胺	4-Chlorophenylethylammonium Iodide	1639014-61-2
4-氰基-2,2,6,6-四甲基哌啶1-氧基 自由基	4-Cyano-2,2,6,6-tetramethylpiperidine 1-Oxyl Free Radical	38078-71-6
4-乙基-4-甲基溴化吗啉	4-Ethyl-4-methylmorpholinium Bromide	65756-41-4
4-氟苯胺氢溴酸盐	4-Fluoroaniline Hydrobromide	85734-18-5
4-氟苯胺盐酸盐	4-Fluoroaniline Hydrochloride	2146/7/8
4-氟苯胺氢碘酸盐	4-Fluoroaniline Hydroiodide	85734-19-6
4-氟苄胺盐酸盐	4-Fluorobenzylamine Hydrochloride	659-41-6
4-氟苯乙胺 盐酸盐	4-Fluorophenethylamine hydrochloride	459-19-8
对氟苯乙胺醋酸	4-Fluorophenylethylammonium Acetate	
4-氟苯甲基溴化胺	4-Fluorophenylmethylammonium Bromide	2270172-94-4
4-氟苯甲基碘化胺	4-Fluorophenylmethylammonium Iodide	2097121-30-5
4-氟苯乙胺三氟醋酸盐	4-Fluorophenylethylammonium Trifluoroacetate	
4-氟苯甲胺三氟醋酸盐	4-Fluorophenylmethylammonium Trifluoroacetate	
4-环氧丙基氧-2,2,6,6-四甲基哌啶1-氧基 自由基	4-Glycidyoxy-2,2,6,6-tetramethylpiperidine 1-Oxyl Free Radical	122413-85-8
4-羟基-2,2,6,6-四甲基哌啶1-氧基苯甲酸酯 自由基 [氧化反应的催化剂]	4-Hydroxy-2,2,6,6-tetramethylpiperidine 1-Oxyl Benzoate Free Radical [Catalyst for Oxidation]	3225-26-1
4-羟基-2,2,6,6-四甲基哌啶1-氧基 自由基	4-Hydroxy-2,2,6,6-tetramethylpiperidine 1-Oxyl Free Radical	2226-96-2
4-羟基苯乙基碘化胺	4-Hydroxy-phenylethylammonium Iodide	1256333-17-1
4-异硫氰酰基-2,2,6,6-四甲基哌啶1-氧基 自由基	4-Isothiocyanato-2,2,6,6-tetramethylpiperidine 1-Oxyl Free Radical	36410-81-8
4-甲基丙烯酰氧-2,2,6,6-四甲基哌啶1-氧自由基	4-Methacryloyloxy-2,2,6,6-tetramethylpiperidine 1-Oxyl Free Radical	15051-46-4
4-甲氧基-2,2,6,6-四甲基哌啶1-氧基 自由基 [氧化反应的催化剂]	4-Methoxy-2,2,6,6-tetramethylpiperidine 1-Oxyl Free Radical [Catalyst for Oxidation]	95407-69-5
4-甲氧基苯乙基溴化胺	4-Methoxyphenylethylammonium Bromide	2705331-53-7
4-甲基-1,3,2-二氧杂硫戊环-2-氧化物 (异构体混合物)	4-Methyl-1,3,2-dioxathiolane 2-Oxide (mixture of isomers)	1469-73-4
4-甲基苯乙胺氢溴酸盐	4-Methyl-2-Phenylethylammonium Bromide	
4-甲基苯乙胺氢碘酸盐	4-Methyl-2-Phenylethylammonium Iodide	
4-甲基苯甲基溴化胺	4-methylphenylmethylammonium bromide	130191-64-9
4-甲基苯甲基碘化胺	4-methylphenylmethylammonium Iodide	289716-73-0
4-氧代-2,2,6,6-四甲基哌啶1-氧 自由基 (升华提纯)	4-Oxo-2,2,6,6-tetramethylpiperidine 1-Oxyl Free Radical (purified by sublimation)	2896-70-0
4-吡啶甲胺氢溴酸盐	4-Pyridinylmethylammonium Bromide	496878-07-0
4-吡啶甲胺氢碘酸盐	4-Pyridinylmethylammonium Iodide	
4-叔丁基苯基溴化胺	4-t-butylphenylammonium Bromide	1446307-72-7
4-叔丁基苯基碘化胺	4-t-butylphenylammonium Iodide	
4-叔丁基苯甲基溴化胺	4-t-butylphenylmethylammonium Bromide	663941-83-1
4-叔丁基苯甲基碘化胺	4-t-butylphenylmethylammonium Iodide	2366873-27-8
4-三氟甲氧基苯乙基碘化胺	4-Trifluoromethoxyphenylethylammonium Iodide	
4-三氟甲基苯胺氢溴酸盐	4-Trifluorophenylammonium Bromide	148819-81-2
4-三氟甲基苯胺氢碘酸盐	4-Trifluorophenylammonium Iodide	
4-三氟甲基苯乙胺氢溴酸盐	4-Trifluorophenylethylammonium Bromide	
4-三氟甲基苯乙胺氢碘酸盐	4-Trifluorophenylethylammonium Iodide	
4-三氟甲基苯甲胺氢溴酸盐	4-Trifluorophenylmethylammonium Bromide	
4-三氟甲基苯甲胺氢碘酸盐	4-Trifluorophenylmethylammonium Iodide	
4-乙烯基-1,3-二氧戊环-2-酮	4-Vinyl-1,3-dioxolan-2-one	4427-96-7

中文名称	英文名称	CAS
5-氨基戊酸氢碘酸盐 (低含水量)	5-Aminovaleric Acid Hydroiodide (Low water content)	1705581-28-7
5-氨基戊酸氢溴酸盐	5-AVABr	2173111-73-2
乙脞氢溴酸盐	Acetamidine Hydrobromide	1040352-82-6
乙脞盐酸盐	Acetamidine hydrochloride	124-42-5
乙脞氢碘酸盐	Acetamidine Hydroiodide	1452099-14-7
己二腈	Adiponitrile	111-69-3
烯丙基甲基砜	Allyl methyl sulfone	16215-14-8
戊基三乙基铵双(三氟甲磺酰)亚胺	Amyltriethylammonium Bis(trifluoromethanesulfonyl)imide	906478-91-9
无水醋酸铅	Anhydrous lead acetate	301-04-2
苯胺氢溴酸盐	Aniline Hydrobromide	542-11-0
苯胺氢碘酸盐	Aniline Hydroiodide	45497-73-2
蒽醌-1,5-二磺酸二钠盐	Anthraquinone-1,5-disulfonic Acid Disodium Salt	853-35-0
苄基(乙基)二甲基铵双(三氟甲磺酰)亚胺	Benzyl(ethyl)dimethylammonium Bis(trifluoromethanesulfonyl)imide	1186103-43-4
苄胺氢溴酸盐	Benzylamine Hydrobromide	37488-40-7
苄胺盐酸盐	Benzylamine Hydrochloride	3287-99-8
苄胺氢碘酸盐(低含水量)	Benzylamine Hydroiodide (Low water content)	45579-91-7
碳酸双(2,2,2-三氟乙基)酯	Bis(2,2,2-trifluoroethyl) Carbonate	1513-87-7
癸二酸双(2,2,6,6-四甲基-4-哌啶基-1-氧基)酯	Bis(2,2,6,6-tetramethyl-4-piperidyl-1-oxyl)	2516-92-9
双(异硫氰酸)(2,2'-二吡啶基-4,4'-二甲酸)(4,4'-二壬基-2,2'-二吡啶基)钌(II)钠盐	Bis(isothiocyanato)(2,2'-bipyridyl-4,4'-dicarboxylato)(4,4'-dinonyl-2,2'-bipyridyl)ruthenium(II) Sodium Salt	871466-65-8
双-PCBM (异构体混合物)	Bis-PCBM (mixture of isomers)	1048679-01-1
丁基醋酸铵	Butanammonium Acetate	1641-36-7
丁胺甲酸	Butanammonium Formate	18721-03-4
丁基胺氢溴酸盐	Butylamine Hydrobromide	15567-09-6
丁胺盐酸盐	Butylamine Hydrochloride	3858-78-4
丁胺氢碘酸盐	Butylamine Hydroiodide	36945-08-1
丁胺铅溴	Butylammonium Lead Bromide	274903-07-0
丁胺铅碘	Butylammonium Lead Iodide	132069-10-4
丁胺甲基磺酸盐	Butylammonium Methyl Sulfonate	53892-58-3
丁胺对甲苯磺酸盐	Butylammonium p-toluenesulfonate	18804-82-5
丁胺三氟醋酸盐	Butylammonium Trifluoroacetate	24922-69-8
丁二胺铅溴	Butyldiammonium Lead Bromide	2712631-79-1
丁二胺铅碘	Butyldiammonium Lead Iodide	1373431-89-0
丁基三甲基铵双(三氟甲磺酰)亚胺	Butyltrimethylammonium Bis(trifluoromethanesulfonyl)imide	258273-75-5
C <sub>60</sub> MC <sub>12</sub>	C <sub>60</sub> MC <sub>12</sub>	403483-19-2
砷化镉	Cadmium arsenide	12006-15-4
双(九氟丁基磺酰基)亚胺钙(II)	Calcium(II) Bis(nonafluorobutanesulfonyl)imide	689282-13-1
碳, 介孔	Carbon, mesoporous	1333-86-4
碳, 介孔	Carbon, mesoporous	1333-86-4
碳, 介孔	Carbon, mesoporous	1333-86-4
碳, 介孔	Carbon, mesoporous	1333-86-4
碳黑	Carbon, mesopor	1333-86-4
碳酸二正丁酯	carbonic acid dibutyl ester	542-52-9
氧化铈, 镨掺杂	Cerium oxide, praseodymium doped	130071-47-5
双(三氟甲磺酰基)亚胺铈(III)	Cerium(III) Bis(trifluoromethanesulfonyl)imide	1046099-39-1
氧化铈(IV) 氧化物	Cerium(IV) oxide	1306-38-3

中文名称	英文名称	CAS
铈(IV)氧化物-钆掺杂	Cerium(IV) oxide-gadolinium doped	
铈(IV)氧化物-钆掺杂	Cerium(IV) oxide-gadolinium doped	
铈(IV)氧化物-钆掺杂	Cerium(IV) oxide-gadolinium doped	
铈(IV)氧化物-钆掺杂	Cerium(IV) oxide-gadolinium doped	
铈(IV)氧化物-钐掺杂	Cerium(IV) oxide-samarium doped	
氧化钐掺杂的氧化铈(IV)	Cerium(IV) oxide-samarium doped	
铈(IV)氧化物-氧化钇掺杂	Cerium(IV) oxide-yttria doped	170033-03-1
标准Z907染料	cis-Bis(isothiocyanato)(2,2'-bipyridyl-4,4'-dicarboxylato)(4,4'-di-nonyl-2'-	502693-09-6
硫化钴(IV)	Cobalt sulfide	12013-10-4
甲基二苯基磷酸酯 (类似物的混合物)	Cresyl Diphenyl Phosphate (mixture of analogue)	26444-49-5
铯铅氯	CsPbCl <sub>3</sub>	15203-83-5
环己胺氢溴酸盐	Cyclohexanemethylamine Hydrobromide	3218-02-8(free)
环己胺氢碘酸盐	Cyclohexanemethylamine Hydroiodide	2153504-15-3
环己胺氢溴酸盐	Cyclohexylamine Hydrobromide	26227-54-3
环己胺氢碘酸盐	Cyclohexylammonium Iodide	45492-87-3
环己基三甲铵双(三氟甲磺酰)亚胺	Cyclohexyltrimethylammonium Bis(trifluoromethanesulfonyl)imide	952155-74-7
环戊胺氢溴酸盐	Cyclopentylammonium Bromide	56383-94-9
环戊胺氢碘酸盐	Cyclopentylammonium Iodide	45372-80-3
环丙胺氢溴酸盐	Cyclopropanammonium Bromide	56383-93-8
环丙胺氢碘酸盐	Cyclopropanammonium Iodide	872-29-7
D 131	D 131	652145-29-4
D 358	D 358	1207638-53-6
癸胺氢溴酸盐	Decaneammonium Bromide	32509-44-7
癸胺甲基磺酸盐	Decaneammonium Methyl Sulfonate	
癸胺对甲苯磺酸盐	Decaneammonium p-toluenesulfonate	126383-73-1
二丁胺氢溴酸盐	Dibutylamine Hydrobromide	10435-44-6
二正丁胺盐酸盐	Dibutylamine hydrochloride	6287-40-7
二丁胺氢碘酸盐	Dibutylamine Hydroiodide	79886-80-9
2,5-二氧杂己二酸二乙酯	Diethyl 2,5-Dioxahexanedioate	35466-87-6
碳酸二乙酯	Diethyl Carbonate	105-58-8
二乙基(2-甲氧基乙基)甲基铵双(氟磺酰)亚胺	Diethyl(2-methoxyethyl)methylammonium Bis(fluorosulfonyl)imide	1079129-48-8
二乙基(甲基)丙基铵双(氟磺酰)亚胺	Diethyl(methyl)propylammonium Bis(fluorosulfonyl)imide	1235234-32-8
二乙胺氢溴酸盐	Diethylamine Hydrobromide	6274-12-0
二乙胺氢碘酸盐	Diethylamine Hydroiodide	19833-78-4
二异丙胺氢碘酸盐	Diisopropylamine Hydroiodide	132396-99-7
2,5-二氧杂己二酸二甲酯	Dimethyl 2,5-Dioxahexanedioate	88754-66-9
碳酸二甲酯(DMC)	Dimethyl carbonate	616-38-6
亚硫酸二甲酯	Dimethyl sulfite	616-42-2
二甲砒	Dimethyl Sulfone	67-71-0
二甲基胺氢溴酸盐	Dimethylamine Hydrobromide	6912-12-5
二甲胺氢碘酸盐	Dimethylamine Hydroiodide	51066-74-1
二甲胺铅碘盐	Dimethylammonium Lead Iodide	1688689-28-2
蒽醌-1,8-二磺酸二钾	Dipotassium Anthraquinone-1,8-disulfonate	14938-42-2
二丙砒	Dipropyl Sulfone	598-03-8
二丙胺氢溴酸盐	Dipropylamine Hydrobromide	7334-96-5
蒽醌-2,6-二磺酸二钠	Disodium Anthraquinone-2,6-disulfonate	853-68-9
DL-1-苯乙基溴化胺	DL-1-Phenylethylammonium Bromide	48104-35-4

中文名称	英文名称	CAS
DL-1-苯乙基碘化胺	DL-1-Phenylethylammonium Iodide	
DNP-Li	DNP-Li	1882849-59-3
十二烷基碘化胺	Dodecylammonium Iodide	34099-97-3
十二胺铅溴	Dodecylammonium Lead Bromide	1239966-53-0
乙醇胺氢溴酸盐	EOABr	23382-12-9
巯基乙胺氢溴酸盐	ESABr	42954-15-4
巯基乙胺氢碘酸盐	ESAI	98650-04-5
盐酸二乙胺	Ethanamine	660-68-4
碳酸甲乙酯	Ethyl methyl carbonate	623-53-0
乙基(2-甲氧基乙基)二甲基铵双(氟磺酰基)亚胺	Ethyl(2-methoxyethyl)dimethylammonium Bis(fluorosulfonyl)imide	1235234-35-1
乙基(2-甲氧基乙基)二甲基铵双(三氟甲磺酰基)亚胺	Ethyl(2-methoxyethyl)dimethylammonium Bis(trifluoromethanesulfonyl)imide	557788-37-1
乙基(3-甲氧基丙基)二甲基铵双(三氟甲磺酰基)亚胺	Ethyl(3-methoxypropyl)dimethylammonium Bis(trifluoromethanesulfonyl)imide	1373334-05-4
乙基(二甲基)(2-苯乙基)铵二(三氟甲磺酰基)亚胺	Ethyl(dimethyl)(2-phenylethyl)ammonium Bis(trifluoromethanesulfonyl)imide	1804970-28-2
乙胺氢溴酸盐	Ethylamine Hydrobromide	593-55-5
乙胺盐酸盐	Ethylamine Hydrochloride	557-66-4
乙胺氢碘酸盐	Ethylamine Hydroiodide	506-58-1
乙胺铅溴	Ethylammonium Lead Bromide	17251-33-1
乙胺铅碘	Ethylammonium Lead Iodide	1640335-04-1
碳酸乙烯酯	Ethylene carbonate	96-49-1
乙二醇双(丙腈)醚	Ethylene Glycol Bis(propionitrile) Ether	3386-87-6
乙二胺二氢溴酸盐	Ethylenediamine Dihydrobromide	624-59-9
乙二胺二氢碘化物	ethylenediamine dihydroiodide	5700-49-2
乙二胺盐酸盐	Ethylenediammonium Dichloride	333-18-6
乙基三甲基碘化铵	Ethyltrimethylammonium Iodide	51-93-4
乙基三苯基碘化膦	Ethyltriphenylphosphonium iodide	4736-60-1
乙基三丙基碘化铵	Ethyltripropylammonium Iodide	15066-80-5
乙基紫精二溴化物	Ethylviologen Dibromide	53721-12-3
氟代碳酸亚乙酯	Fluoroethylene carbonate	114435-02-8
甲脒甲酸	Formamidine formate	116013-65-1
甲脒氢溴酸盐(低含水量)	Formamidine Hydrobromide (Low water content)	146958-06-7
甲脒盐酸盐	Formamidine hydrochloride	6313-33-3
甲脒氢碘酸盐(低含水量)	Formamidine Hydroiodide (Low water content)	879643-71-7
甲脒硫氰酸盐	Formamidine Thiocyanate	1821033-48-0
甲脒铅溴盐	Formamidinium Lead Bromide	1008105-17-6
甲脒铅溴碘盐	Formamidinium Lead Bromide Iodide	1616115-27-5
甲脒铅氯碘盐	Formamidinium Lead Chloride Iodide	1616115-25-3
甲脒铅碘盐	Formamidinium Lead Iodide	1451592-07-6
四氟硼酸甲脒	Formamidinium Tetrafluoroborate	
反丁烯二腈	Fumaronitrile	764-42-1
加尔万氧基自由基	Galvinoxyl Free Radical	2370-18-5
戊二腈	Glutaronitrile	544-13-8
丙三醇1,2-碳酸酯	Glycerol 1,2-Carbonate	931-40-8
石墨烯管	Graphene tube	
胍基氢溴酸盐	Guanidinium Bromide	19244-98-5
胍基氢碘酸盐	Guanidinium Iodide	19227-70-4
己基醋酸铵	HAAC	23239-72-7
己胺铅溴	Hexylammonium Lead Bromide	274903-09-2



中文名称	英文名称	CAS
六氟磷酸锂溶液	Lithium hexafluorophosphate solution	
六氟磷酸锂溶液	Lithium hexafluorophosphate solution	
六氟磷酸锂溶液	Lithium hexafluorophosphate solution	
六氟磷酸锂溶液	Lithium hexafluorophosphate solution	
六氟磷酸锂溶液	Lithium hexafluorophosphate solution	
碘化锂 (99.9%, 微量金属品) (低含水量)	Lithium Iodide (99.9%, trace metals basis) (Low water content)	10377-51-2
磷酸铁锂	Lithium iron(II) phosphate	15365-14-7
二氧化锰锂	Lithium manganese dioxide	12162-79-7
锂锰镍氧化物	Lithium manganese nickel oxide	
锂锰镍氧化物	Lithium manganese nickel oxide	12031-75-3
锂锰氧化物	Lithium manganese oxide	
锂锰氧化物	Lithium manganese oxide	12057-17-9
磷酸锰(II)锂	Lithium Manganese(II) Phosphate	13826-59-0
偏硼酸锂	Lithium Metaborate	13453-69-5
偏磷酸锂	Lithium Metaphosphate	13762-75-9
硅酸锂	Lithium Metasilicate	10102-24-6
锂镍钴铝氧化物	Lithium nickel cobalt aluminium oxide	
锂镍钴铝氧化物	Lithium nickel cobalt aluminium oxide	193214-24-3
锂镍钴氧化物	Lithium nickel cobalt oxide	113066-89-0
二氧化镍锂	Lithium nickel dioxide	12031-65-1
锂镍锰钴氧化物	Lithium nickel manganese cobalt oxide	
锂镍锰钴氧化物	Lithium nickel manganese cobalt oxide	346417-97-8
锂镍锰钴氧化物 (LiNi <sub>0.4</sub> Mn <sub>0.4</sub> Co <sub>0.2</sub> O <sub>2</sub> )	Lithium Nickel Manganese Cobalt Oxide (LiNi <sub>0.4</sub> Mn <sub>0.4</sub> Co <sub>0.2</sub> O <sub>2</sub> )	193215-96-2
锂镍锰钴氧化物 (LiNi <sub>0.6</sub> Mn <sub>0.2</sub> Co <sub>0.2</sub> O <sub>2</sub> )	Lithium Nickel Manganese Cobalt Oxide (LiNi <sub>0.6</sub> Mn <sub>0.2</sub> Co <sub>0.2</sub> O <sub>2</sub> )	193215-05-3
锂镍锰钴氧化物 (LiNi <sub>0.8</sub> Mn <sub>0.1</sub> Co <sub>0.1</sub> O <sub>2</sub> )	Lithium Nickel Manganese Cobalt Oxide (LiNi <sub>0.8</sub> Mn <sub>0.1</sub> Co <sub>0.1</sub> O <sub>2</sub> )	179802-95-0
磷酸锂镍(II)	Lithium Nickel(II) Phosphate	13977-83-8
铌酸锂 (99.99%, 微量金属品)	Lithium Niobate (99.99%, trace metals basis)	12031-63-9
九氟-1-丁烷磺酸锂	Lithium Nonafluoro-1-butanefluorobutanesulfonate	131651-65-5
氧化锂	Lithium Oxide	12057-24-8
磷酸锂	Lithium Phosphate	10377-52-3
二氟磷酸锂	Lithium Phosphorodifluoridate	24389-25-1
丙酮酸锂 一水合物 [GOT, GPT标准品的保护剂]	Lithium Pyruvate Monohydrate [Guaranteed for Standard to GOT, GPT]	2922-61-4
硬脂酸锂	Lithium Stearate	4485-12-5
硫酸锂	Lithium Sulfate	10377-48-7
硫化锂	Lithium Sulfide	12136-58-2
钽酸锂 (99.99%, 微量金属品)	Lithium Tantalate (99.99%, trace metals basis)	12031-66-2
四硼酸锂	Lithium Tetraborate	12007-60-2
四氟硼酸锂	Lithium tetrafluoroborate	14283-07-9
钛酸锂	Lithium titanate	
钛酸锂/锂钛氧化物 (尖晶石)	Lithium Titanate (Spinel)	12031-95-7
磷酸钛锂 (LiTi <sub>2</sub> (PO <sub>4</sub> ) <sub>3</sub> )	Lithium Titanium Phosphate (LiTi <sub>2</sub> (PO <sub>4</sub> ) <sub>3</sub> )	30622-39-0
三氰基甲烷化锂-乙腈络合物	Lithium Tricyanomethanide - Acetonitrile Complex	
三氰基甲烷化锂-乙腈络合物	Lithium Tricyanomethanide - Acetonitrile Complex	
三氟甲磺酸锂	Lithium Trifluoromethanesulfonate	33454-82-9
磷酸钒锂 (Li <sub>3</sub> V <sub>2</sub> (PO <sub>4</sub> ) <sub>3</sub> )	Lithium Vanadium Phosphate (Li <sub>3</sub> V <sub>2</sub> (PO <sub>4</sub> ) <sub>3</sub> )	84159-18-2
甲基甲酸铵	MAFA	25596-28-5

中文名称	英文名称	CAS
丙二腈	Malonitrile	109-77-3
丙二腈标准溶液	Malonitrile	109-77-3(methanol)
二氧化锰	Manganese dioxide	1313-13-9
甲基丙酸铵	MAPA	54856-22-3
熔体流动速率标准物质(聚丙烯)	Melt Flow Rate (PP)	9003-07-0
MeO-2PACz	MeO-2PACz	2377770-18-6
甲氧基乙腈	Methoxyacetonitrile	1738-36-9
2,2,2-三氟乙基碳酸甲酯	Methyl 2,2,2-Trifluoroethyl Carbonate	156783-95-8
甲胺氰酸盐	Methylamine Cyanate	63405-91-4
甲胺氢溴酸盐	Methylamine Hydrobromide	6876-37-5
甲胺盐酸盐	Methylamine hydrochloride	593-51-1
甲基溴化胺	Methylammonium Bromide	6816-37-5
甲铵六氟磷酸盐	Methylammonium Hexafluorophosphate	28302-50-3
甲铵铅溴碘盐	Methylammonium Lead Bromide Iodide	1472068-58-8
甲铵铅溴碘盐	Methylammonium Lead Bromide Iodide	1446121-02-3
甲铵铅氯溴盐	Methylammonium Lead Chloride Bromide	479066-07-4
甲铵铅氯溴盐	Methylammonium Lead Chloride Bromide	479066-09-6
甲铵铅氯碘盐	Methylammonium Lead Chloride Iodide	1446121-07-8
甲铵铅氯碘盐	Methylammonium Lead Chloride Iodide	1472068-56-6
辛胺甲基磺酸盐	Methylammonium Methyl Sulfonate	1450981-06-2
甲铵四氟硼酸盐	Methylammonium Tetrafluoroborate	42539-74-2
硫氰酸甲铵	Methylammonium Thiocyanate	61540-63-4
甲铵三溴铅酸盐	Methylammonium tribromoplumbate	69276-13-7
亚甲二胺二氢溴酸盐	Methyldiammonium dibromide	
亚甲二胺二盐酸盐	Methyldiammonium dichloride	57166-92-4
亚甲二胺二碘酸盐	Methyldiammonium diiodide	
甲基三正辛铵合双(三氟甲烷磺酰)亚胺	Methyltri-n-octylammonium Bis(trifluoromethanesulfonyl)imide	375395-33-8
吗啉氢溴酸盐	Morpholine Hydrobromide	6377-82-8
吗啉氢碘酸盐	Morpholine Hydroiodide	58464-45-2
N,N'-双(三-p-甲苯基磷酰亚烷基)苯-1,4-二胺	N,N'-Bis(tri-p-tolylphosphoranylidene)benzene-1,4-diamine	51870-56-5
N,N-二乙基-1,2-乙二胺二盐酸盐	N,N-Diethyl-1,2-ethanediamine Dihydrochloride	52198-62-6
N,N-二乙基-N-甲氧基-N-甲基乙烷-1-铵双[(三氟甲基)磺酰基]酰胺	N,N-Diethyl-2-methoxy-N-methylethan-1-aminium Bis[(trifluoromethyl)sulfonyl]amide	464927-84-2
N,N-二甲基-1,3-丙二胺二氢溴酸盐	N,N-Dimethyl-1,3-propanediamine	
N,N-二甲基-1,3-丙二胺氢溴酸盐	N,N-Dimethyl-1,3-propanediammonium Dibromide	2710685-13-3
N,N-二甲基-1,3-丙二胺盐酸盐	N,N-Dimethyl-1,3-propanediammonium Dichloride	
N,N-二甲基-1,3-丙二胺氢碘酸盐	N,N-Dimethyl-1,3-propanediammonium Diiodide	2561497-43-4
N,N-二甲基-1,3-丙二胺铅溴	N,N-Dimethyl-1,3-propanediammonium Lead	
N,N-二甲基-1,3-丙二胺铅碘	N,N-Dimethyl-1,3-propanediammonium Lead	
N,N-二甲基乙二胺二氢溴酸盐	N,N-Dimethylethylenediamine Dihydrobromide	1245570-04-0
N,N-二甲基乙二胺二氢碘酸盐	N,N-Dimethylethylenediamine Dihydroiodide	244234-52-4
N2,N2',N7,N7'-四(4-甲氧基苯基)-N2,N2',N7,N7'-四(2-(9,9-二甲基芴))螺[芴-9,9'-氧杂蒽]-2,2',7,7'-四胺	N7,N7'-trans(4-methoxyphenyl)-N2,N2',N7,N7'-trans(9,9-dimethylfluorene-2-yl)spiro[fluorene-9,9'-xanthene]-2,2',7,7'-tetraamine	2128308-79-0
N749黑染料	N749 Black Dye	359415-47-7
Nafion™ 全氟化树脂溶液(D2020CS)	Nafion™ perfluorinated resin solution	31175-20-9(Water)
新戊胺氢溴酸盐	Neopentylamine Hydrobromide	

中文名称	英文名称	CAS
新戊胺氢碘酸盐	Neopentylamine Hydroiodide	
氧化镍 - 钪稳定氧化锆, NiO/ScSZ-10	Nickel oxide - scandia-stabilized zirconia,	
氧化镍 - 钪稳定氧化锆, NiO/ScSZ-6	Nickel oxide - scandia-stabilized zirconia,	
氧化镍 - 氧化钇稳定的氧化锆	Nickel oxide - Ytria-stabilized zirconia	
过氧化镍	Nickel(II) peroxide	86676-91-7
氮掺杂有序介孔碳	Nitrogen doped ordered mesoporous carbon	
正辛胺氢溴酸盐	n-Octylamine Hydrobromide	14846-47-0
正辛胺盐酸盐	n-Octylamine Hydrochloride	142-95-0
壬基溴化铵	Nonylammonium Bromide	43146-93-6
壬基碘化铵	Nonylammonium Iodide	63478-97-7
辛基醋酸铵	OAAc	2016-40-2
十八胺氢溴酸盐	Octadecanammium Bromide	68810-17-3
十八胺氢碘酸盐	Octadecanammium iodide	17535-32-9
辛胺铅溴	Octylammium Lead Bromide	882567-53-5
辛胺铅碘	Octylammium Lead Iodide	
辛胺对甲苯磺酸盐	Octylammium p-toluenesulfonate	61165-94-4
9-十八烯基溴化胺 (油胺溴)	Oleylammonium Bromide OAmBr	2024602-96-6
9-十八烯基碘化胺 (油胺碘)	Oleylammonium Iodide OAmI	1802520-56-4
有序介孔碳	Ordered mesoporous carbons	
苯乙基醋酸铵	PEAAc	24722-38-1
戊基醋酸铵	Pentylammium Acetate	64370-74-7
苯丁基醋酸铵	PhBAAc	24722-33-6
苯乙基氯化胺	Phenethylammium Chloride	156-28-5
苯丁基氯化胺	Phenylbutanammium Chloride	30684-06-1
苯丁基碘化胺	Phenylbutanammium iodide	
苯乙胺铅溴	Phenylethylammium Lead Bromide	131457-16-4
苯乙胺铅碘	Phenylethylammium Lead Iodide	131457-08-4
苯乙胺三氟醋酸盐	Phenylethylammium Trifluoroacetate	203525-17-1
苯甲胺铅溴	Phenylmethylammium Lead Bromide	17251-04-6
苯甲胺铅碘	Phenylmethylammium Lead Iodide	256222-08-9
苯丙基碘化胺	Phenylpropylammium Iodide	2673196-34-2
物理法石墨烯粉末	Physical method of graphene powder	
物理法石墨烯粉末	Physical method of graphene powder	
物理法石墨烯粉末	Physical method of graphene powder	
物理法石墨烯粉末	Physical method of graphene powder	
物理法石墨烯粉末	Physical method of graphene powder	
哌嗪二氢碘酸盐	Piperazine Dihydriodide	58464-47-4
哌嗪二氢溴酸盐	Piperazine Dihydrobromide	59813-05-7
哌嗪单溴	Piperazinehydrobromide	59813-07-9
哌嗪单碘	Piperazinehydroiodide	56310-12-4
聚(9,9-n-二辛基-2,7-芴-alt-9-异辛基-3,6-咔唑)	Poly(9,9-n-dioctyl-2,7-fluorene-alt-9-iso-octyl-3,6-carbazole)	
聚四氟乙烯	Poly(tetrafluoroethylene)	9002-84-0
聚四氟乙烯	Poly(tetrafluoroethylene)	9002-84-0
聚四氟乙烯微粉	Poly(tetrafluoroethylene)	9002-84-0
聚氯乙烯	Poly(vinyl chloride)	9002-86-2
聚氯乙烯	Poly(vinyl chloride)	9002-86-2
聚氯乙烯	Poly(vinyl chloride)	9002-86-2
聚氯乙烯	Poly(vinyl chloride)	9002-86-2
聚氯乙烯	Poly(vinyl chloride)	9002-86-2
聚[[4,4'-(N-(4-仲丁基苯基)二苯胺)]	Poly[4,4'-(N-(4-secbutylphenyl)diphenylamine]	860495-18-7



中文名称	英文名称	CAS
聚[双(4-苯基)(2,4-二甲基苯基)胺]	Poly[bis(4-phenyl)(2,4-dimethylphenyl)amine]	313996-10-0
聚丙烯	Polypropylene	9003-07-0
聚丙烯	Polypropylene	9003-07-0
聚丙烯	Polypropylene	9003-07-0
聚丙烯	Polypropylene	9003-07-0
聚丙烯	Polypropylene	9003-07-0
聚丙烯	Polypropylene	9003-07-0
聚丙烯	Polypropylene	9003-07-0
聚丙烯	Polypropylene	9003-07-0
聚四氟乙烯微粉	Polytetrafluoroethylene	9002-84-0( powder)
(氟磺酰)(三氟甲磺酰)亚氨基钾	Potassium (Fluorosulfonyl)(trifluoromethanesulfonyl)imide	860653-59-4
双(三氟甲基磺酰基)酰亚胺钾	Potassium Bis(trifluoromethanesulfonyl)imide	90076-67-8
双九氟-1-丁烷磺酰亚胺钾	Potassium Bisnonafluoro-1-butanesulfonimidate	129135-87-1
六氟磷酸钾	Potassium Hexafluorophosphate	17084-13-8
四氟硼酸钾	Potassium Tetrafluoroborate	14075-53-7
三氰基甲烷化钾	Potassium Tricyanomethanide	34171-69-2
三(三氟甲烷磺酰基)甲基化钾	Potassium Tris(trifluoromethanesulfonyl)methanide	114395-69-6
聚丙烯熔体流动速率标准物质	Ppolypropylene(melt flow rate)	9003-07-0
丙胺氢溴酸盐	Propylamine Hydrobromide	4905-83-3
丙胺盐酸盐	Propylamine Hydrochloride	556-53-6
丙胺氢碘酸盐	Propylamine Hydroiodide	14488-45-0
丙胺铅溴	Propylammonium Lead Bromide	2338841-75-9
丙胺铅碘	Propylammonium Lead Iodide	1646171-78-9
碳酸丙烯酯	Propylene carbonate	108-32-7
吡啶氢碘酸盐	Pyridine Hydroiodide	18820-83-2
吡咯烷氢碘酸盐	Pyrrolidine Hydriodide	45361-12-4
r-1-甲基萘甲胺溴	r-(-)-1-(1-Naphthyl)ethylammonium Bromide	2484838-97-1
r-1-甲基萘甲胺碘	r-(+)-1-(1-Naphthyl)ethylammonium Iodide	
r-甲基苄胺溴	r-(+)-1-Phenylethylammonium Bromide	48104-28-5
r-甲基苄胺碘	r-(+)-1-PhenylethylammoniumIodide	
r-p-Br-甲基苄胺碘	r-(+)-4-Br- $\alpha$ -MethylbenzylammoniumIodide	
s-1-甲基萘甲胺溴	s-(-)-1-(1-Naphthyl)ethylammonium Bromide	2484838-96-0
s-1-甲基萘甲胺碘	s-(-)-1-(1-Naphthyl)ethylammonium Iodide	
s-(-)-2-甲基萘甲基溴化胺	s-(-)-2-(1-Naphthyl)ethylammonium Bromide	
s-(-)-2-甲基萘甲基碘化胺	s-(-)-2-(1-Naphthyl)ethylammonium Iodide	
s-甲基苄胺溴	s-(-)-1-Phenylethylammonium Bromide	48104-32-1
s-甲基苄胺碘	s-(-)-1-PhenylethylammoniumIodide	
s-p-Br-甲基苄胺溴	s-(-)-4-Br- $\alpha$ -Methylbenzylammonium Bromide	
s-p-Br-甲基苄胺碘	s-(-)-4-Br- $\alpha$ -Methylbenzylammonium Iodide	
蒽醌-1-磺酸钠	Sodium Anthraquinone-1-sulfonate	128-56-3
双(氟磺酰)亚胺钠	Sodium Bis(fluorosulfonyl)imide	100669-96-3
双(九氟丁基磺酰基)亚胺钠	Sodium Bis(nonafluorobutanesulfonyl)imide	129135-86-0
氟硼酸钠	Sodium fluoroborate	13755-29-8
六氟磷酸钠	Sodium Hexafluorophosphate	21324-39-0
三氰基甲烷化钠	Sodium Tricyanomethanide	36603-80-2
丁二腈	Succinonitrile	110-61-2
环丁砜	Sulfolane	126-33-0
叔丁基胺氢溴酸盐	tert-Butylamine Hydrobromide	60469-70-7
叔丁基胺氢碘酸盐	tert-Butylamine Hydroiodide	39557-45-4

中文名称	英文名称	CAS
四戊基溴化铵	Tetraamylammonium Bromide	866-97-7
四戊基氯化铵	Tetraamylammonium Chloride	4965-17-7
四戊基碘化铵	Tetraamylammonium Iodide	2498-20-6
四丁基乙酸铵	Tetrabutylammonium Acetate	10534-59-5
四丁基溴化铵 (TBAB)	Tetrabutylammonium Bromide	1643-19-2
四丁基氯化铵	Tetrabutylammonium Chloride	1112-67-0
四丁基六氟磷酸铵	Tetrabutylammonium Hexafluorophosphate	3109-63-5
四丁基碘化铵	Tetrabutylammonium Iodide	311-28-4
对甲苯磺酸四丁铵	Tetrabutylammonium p-Toluenesulfonate	7182-86-7
四丁基四氟硼酸铵	Tetrabutylammonium Tetrafluoroborate	429-42-5
四丁基三氟甲磺酸铵	Tetrabutylammonium Trifluoromethanesulfonate	35895-70-6
四丁基溴化磷	Tetrabutylphosphonium Bromide	3115-68-2
四丁基六氟磷酸磷	Tetrabutylphosphonium Hexafluorophosphate	111928-21-3
四丁基 O,O-二乙基二硫代磷酸磷	Tetrabutylphosphonium O,O-Diethyl Phosphorodithioate	96131-57-6
四丁基四氟硼酸磷	Tetrabutylphosphonium Tetrafluoroborate	1813-60-1
四氰基乙烯 (升华提纯)	Tetracyanoethylene (purified by sublimation)	670-54-2
四乙基碘化铵	Tetraethylammonium iodide	68-05-3
四庚基溴化铵	Tetraheptylammonium Bromide	4368-51-8
四庚基碘化铵	Tetraheptylammonium Iodide	3535-83-9
四己基溴化铵	Tetrahexylammonium Bromide	4328-13-6
四己基碘化铵	Tetrahexylammonium Iodide	2138-24-1
四甲基碘化铵	Tetramethylammonium iodide	75-58-1
四正辛基溴化铵	Tetra-n-octylammonium Bromide	14866-33-2
四正辛基碘化铵	Tetra-n-octylammonium Iodide	16829-91-7
四正辛基溴化磷	Tetra-n-octylphosphonium Bromide	23906-97-0
四苯基碘化磷	Tetraphenylphosphonium Iodide	2065-67-0
四丙基氯化铵	Tetrapropylammonium Chloride	5810-42-4
四丙基碘化铵	Tetrapropylammonium iodide	631-40-3
溴化锡(II)	Tin(II) bromide	10031-24-0
磷酸三烯丙酯	Triallyl Phosphate	1623-19-4
磷酸三戊酯	Triamyl Phosphate	2528-38-3
三丁基(2-甲氧基乙基)磷双(三氟甲磺酰基)亚胺	Tributyl(2-methoxyethyl)phosphonium Bis(trifluoromethanesulfonyl)imide	959698-44-3
三丁基(乙基)磷二乙基磷酸盐	Tributyl(ethyl)phosphonium Diethyl Phosphate	20445-94-7
三丁基(甲基)铵二氰胺盐	Tributyl(methyl)ammonium Dicyanamide	1262230-03-4
三丁基(甲基)磷磷酸二甲酯盐	Tributyl(methyl)phosphonium Dimethyl Phosphate	20445-88-9
三丁基十六烷基溴化磷	Tributylhexadecylphosphonium Bromide	14937-45-2
三丁基己基溴化磷	Tributylhexylphosphonium Bromide	105890-71-9
三丁基甲基铵双(三氟甲基磺酰基)亚胺	Tributylmethylammonium bis(trifluoromethylsulfonyl)imide	405514-94-5
三丁基甲基碘化铵	Tributylmethylammonium Iodide	3085-79-8
三丁基甲磷双(三氟甲磺酰基)亚胺	Tributylmethylphosphonium Bis(trifluoromethanesulfonyl)imide	324575-10-2
三丁基甲基碘化磷	Tributylmethylphosphonium Iodide	1702-42-7
三丁基正辛基溴化磷	Tributyl-n-octylphosphonium Bromide	57702-65-5
三丁基碘化铊	Tributylsulfonium Iodide	18146-62-8
磷酸三甲苯酯 (异构体混合物)	Tricresyl Phosphate (mixture of isomers)	1330-78-5
磷酸三乙酯	Triethyl Phosphate	78-40-0
三乙基苯基碘化铵	Triethylphenylammonium Iodide	1010-19-1
三乙基铊双(三氟甲基磺酰基)亚胺	Triethylsulfonium	321746-49-0

中文名称	英文名称	CAS
硼酸三己酯	Trihexyl Borate	5337-36-0
三己基(十四烷基)磷双((三氟甲基)磺酰基)酰胺	Trihexyl(tetradecyl)phosphonium Bis((trifluoromethyl)sulfonyl)amide	460092-03-9
三己基(十四烷基)氯化磷	Trihexyl(tetradecyl)phosphonium Chloride	258864-54-9
三己基(十四烷基)磷二氰胺盐	Trihexyl(tetradecyl)phosphonium Dicyanamide	701921-71-3
磷酸三间甲苯酯	Tri-m-cresyl Phosphate	563-04-2
硼酸三甲酯	Trimethyl Borate	121-43-7
磷酸三甲酯	Trimethyl Phosphate	512-56-1
三甲基苯基碘化铵	Trimethylphenylammonium Iodide	98-04-4
三甲基丙基铵双(三氟甲烷磺酰)亚胺	Trimethylpropylammonium Bis(trifluoromethanesulfonyl)imide	268536-05-6
三甲基碘化铊	Trimethylsulfonium Iodide	2181-42-2
磷酸三邻甲苯酯	Tri-o-cresyl Phosphate	78-30-8
三对甲苯基磷酸酯	Tri-p-cresyl Phosphate	78-32-0
磷酸三苯酯	Triphenyl Phosphate	115-86-6
磷酸三(1,1,1,3,3,3-六氟-2-丙基)酯	Tris(1,1,1,3,3,3-hexafluoro-2-propyl) Phosphate	66489-68-7
三(1,10-菲咯啉)双(六氟磷酸盐)钴(II)	Tris(1,10-phenanthroline)cobalt(II) bis(hexafluorophosphate)	31876-74-1
三(1,10-菲咯啉)三(六氟磷酸盐)钴(III)	Tris(1,10-phenanthroline)cobalt(III) tris(hexafluorophosphate)	28277-59-0
磷酸三(1,3-二氯-2-丙基)酯	Tris(1,3-dichloro-2-propyl) Phosphate	13674-87-8
磷酸三(1H,1H,5H-八氟戊基)酯	Tris(1H,1H,5H-octafluoropentyl) Phosphate	355-86-2
硼酸三(2,2,2-三氟乙基)酯	Tris(2,2,2-trifluoroethyl) Borate	659-18-7
磷酸三(2,2,2-三氟乙基)酯	Tris(2,2,2-trifluoroethyl) Phosphate	358-63-4
亚磷酸三(2,2,2-三氟乙基)酯	Tris(2,2,2-trifluoroethyl) Phosphite	370-69-4
磷酸三(2-正丁氧乙基)酯	Tris(2-butoxyethyl) Phosphate	78-51-3
磷酸三(2-氯乙基)酯	Tris(2-chloroethyl) Phosphate	115-96-8
磷酸三(2-乙基己基)酯	Tris(2-ethylhexyl) Phosphate	78-42-2
三[4-(5-苯基噻吩-2-基)苯基]胺	Tris[4-(5-phenylthiophen-2-yl)phenyl]amine	803727-09-5
超高电容量多孔碳活性炭	Ultrahigh capacitance porous carbon activated	
超高电容量多孔碳活性炭	Ultrahigh capacitance porous carbon activated	
V886	V886	1801701-58-5
碳酸亚乙烯酯(含稳定剂BHT)	Vinylene Carbonate (stabilized with BHT)	872-36-6
氯化锌(高纯度 & 低含水量)	Zinc Chloride (High Purity & Low water content)	7646-85-7
三氟甲烷磺酸锌(II)[电解质用]	Zinc(II) Trifluoromethanesulfonate [for	54010-75-2
氧化锆(IV)氧化钪稳定	Zirconium(IV) oxide-scandia stabilized	
氧化锆(IV)氧化钪稳定	Zirconium(IV) oxide-scandia stabilized	151575-30-3
氧化锆(IV)-氧化钇稳定	Zirconium(IV) oxide-yttria stabilized	
氧化锆(IV)-氧化钇稳定	Zirconium(IV) oxide-yttria stabilized	
氧化锆(IV)-氧化钇稳定	Zirconium(IV) oxide-yttria stabilized	113482-02-3
氧化锆(IV)-氧化钇稳定	Zirconium(IV) oxide-yttria stabilized	114168-16-0
α-八噻吩	α-Octithiophene	113728-71-5
α-四联噻吩	α-Quaterthiophene	5632-29-1
α-喹噻吩	α-Quinquethiophene	5660-45-7
α-七噻吩	α-Septithiophene	86100-63-2
α-西噻吩	α-Sexithiophene	88493-55-4
碘化钾胺	Methylammonium Iodide	14965-49-2
环丁烷四甲酸二酐	Cyclobutane-1,2,3,4-tetracarboxylic dianhydride	
丙氧化新戊二醇二丙烯酸酯		84170-74-1
二丙二醇二丙烯酸酯		57472-68-1
二缩三丙二醇二丙烯酸酯		42978-66-5

中文名称	英文名称	CAS
1,6-己二醇二丙烯酸酯		13048-33-4
丙烯酸吗啉		5117-12-4
三羟甲基丙烷三丙烯酸酯		15625-89-5
新戊二醇二丙烯酸酯		2223-82-7
乙氧基化三羟甲基丙烷三丙烯酸酯		
六甲基二硅氮烷	Hexamethyldisilazane	999-97-3
1-羟基环己基苯基甲酮	1-Hydroxycyclohexyl phenyl ketone	947-19-3
安息香二甲醚	2,2-Dimethoxy-2-phenylacetophenone	24650-42-8
2-羟基-2-甲基-1-苯基-1-丙酮	2-Hydroxy-2-methylpropiophenone	7473-98-5
2-异丙基硫杂蒽酮 (2、4-异构体混合物)	2-Isopropylthioxanthone	5495-84-1
4-氯二苯甲酮	4-Chlorobenzophenone	134-85-0
邻苯甲酰苯甲酸甲酯	Methyl 2-benzoylbenzoate	606-28-0
2,4,6-(三甲基苯甲酰基)二苯基氧化膦		75980-60-8
2,4,6-三甲基苯甲酰基苯基膦酸乙酯		84434-11-7
2-甲基-1-(4-甲硫基苯基)-2-吗啉基-1-丙酮		71868-10-5
4-苯基二苯甲酮	4-Benzoylbiphenyl	2128-93-0
4-二甲氨基苯甲酸乙酯	Ethyl 4-dimethylaminobenzoate	10287-53-3
氨增效剂		21245-02-3
二苯甲酮及衍生物		90-93-7
二苯甲酮及衍生物		119-61-9
二苯甲酮及衍生物		134-84-9
可见光PI		125051-32-3
硫杂蒽酮类		82799-44-8
$\alpha$ -羟基酮		106797-53-9
$\alpha$ -羟基酮		119313-12-1
$\alpha$ -羟基酮		119344-86-4
苯甲酰甲酸酯类		15206-55-0
酰基磷氧化物		162881-26-7
试剂Cerium(IV) oxide-samarium doped	Cerium(IV) oxide-samarium doped	
试剂Lanthanum nickelate, strontium doped	Lanthanum nickelate, strontium doped	
试剂Lanthanum strontium manganite	Lanthanum strontium manganite	
试剂Zirconium(IV) oxide-scandia stabilized	Zirconium(IV) oxide-scandia stabilized	
试剂Zirconium(IV) oxide-yttria stabilized	Zirconium(IV) oxide-yttria stabilized	
乙二醇二甲醚	1,2-Dimethoxyethane	110-71-4
甲烷二磺酸亚甲酯	1,5,2,4-Dioxadithiane 2,2,4,4-Tetraoxide	99591-74-9
二乙二醇二甲醚	2-Methoxyethyl Ether	111-96-6
双(三氟甲磺酰基)亚胺钡(II)	Barium(II) Bis(Trifluoromethanesulfonyl)Imide	168106-22-7
双(三氟甲磺酰基)亚胺锌(II)	Bis(Trifluoromethanesulfonyl)Imide Zinc(II) Salt	168106-25-0
双(氟磺酰)亚胺钙	Calcium Bis(Fluorosulfonyl)Imide	1906900-37-5
双(三氟甲磺酰基)酰亚胺钙(II)	Calcium(II) Bis(Trifluoromethanesulfonyl)Imide	165324-09-4
钷掺杂的氧化铈(IV)	Cerium(IV) oxide-gadolinium doped	
环己基苯	Cyclohexylbenzene	827-52-1
石墨	Graphite	7782-42-5
石墨粉	Graphite	7782-42-5,(稳定提供多瓶与
钴铁酸锶镧	Lanthanum strontium cobalt ferrite	
锰酸锶镧	Lanthanum strontium manganite	
(氟磺酰)(三氟甲磺酰)亚胺基锂	Lithium (Fluorosulfonyl)(Trifluoromethanesulfonyl)Imide	192998-62-2
乙酸锂二水合物	Lithium Acetate Dihydrate	6108-17-4

中文名称	英文名称	CAS
氧化铝锂	Lithium Aluminum Oxide	12003-67-7,(稳定提供多瓶与大包装)
砷酸锂	Lithium Arsenate	13478-14-3
双(五氟乙基磺酰基)亚氨基锂	Lithium Bis(Pentafluoroethanesulfonyl)Imide	132843-44-8
双三氟甲烷磺酰亚胺锂	Lithium bis(trifluoromethanesulphonyl)imide	90076-65-6
溴化锂二水合物	Lithium Bromide Hydrate	13453-70-8
溴化锂水合物	Lithium Bromide Monohydrate	85017-82-9
氯化锂, 无水	Lithium Chloride	7447-41-8
氯化锂单水合物	Lithium Chloride Hydrate	16712-20-2
氯化锂x水合物	Lithium Chloride Xhydrate	85144-11-2
氯化钾锂	Lithium Chloride-Potassium Chloride	65567-96-6
磷酸二氢锂	Lithium Dihydrogenphosphate	13453-80-0
氟化锂	Lithium Fluoride	7789-24-4
甲酸锂一水合物	Lithium Formate Monohydrate	6108-23-2
氢氧化锂, 无水	Lithium Hydroxide	1310-65-2
氢氧化锂一水合物	Lithium hydroxide monohydrate	1310-66-3
碘酸锂	Lithium Iodate	13765-03-2
碘化锂水合物	Lithium Iodide Hydrate	85017-80-7
铁酸锂	LITHIUM IRON(III) OXIDE 95	12022-46-7
偏硼酸锂二水合物	Lithium Metaborate Dihydrate	15293-74-0
偏镓酸锂	Lithium metagallate	944251-30-3
钼酸锂	Lithium Molybdate	13568-40-6
硝酸锂	Lithium Nitrate	7790-69-4
草酸锂	Lithium Oxalate	553-91-3
高碘酸锂二水合物, Reagent Grade	Lithium Periodate Dihydrate, Reagent Grade	21111-84-2
四氯金(III)酸锂	Lithium Tetrachloroaurate(III)	3145-91-3
四氯钋(II)酸锂	Lithium Tetrachloropalladate(II)	15525-45-8
硫氰酸锂水合物	Lithium Thiocyanate Hydrate	123333-85-7
钛酸锂	Lithium titanium oxide	12031-82-2
钒酸锂 (metals basis)	Lithium Vanadium Oxide (Metals Basis)	15060-59-0
锆酸锂	Lithium Zirconate	12031-83-3
钨酸锂	Lithiumtungstate	13568-45-1
双(三氟甲磺酰基)酰亚胺镁(II)	Magnesium Bis(Trifluoromethane Sulfone)Imide	133395-16-1
双(三氟甲磺酰基)亚胺锰(II)	Manganese(Ii) Bis(Trifluoromethanesulfonyl)Imide	207861-55-0
N-甲基吡咯烷酮	N-Methyl-2-Pyrrolidinone	872-50-4
双氟磺酰亚胺钾盐	Potassium Bis(Fluorosulfonyl)Imide	14984-76-0
双(三氟甲基磺酰基)酰亚胺钠	Sodium Bis(Trifluoromethanesulfonyl)Imide	91742-21-1
三氟甲磺酸钠	Sodium Triflate	2926-30-9
六氟异丙氧基环三磷腈	2,2,4,4,6,6-hexahydro-2,2,4,4,6,6-hexakis[2,2,2-trifluoro-1-(trifluoromethyl)ethoxy]-1,3,5,2,4,6-	80192-24-1
2,2,4,4,6-五氟-6-苯氧基-2λ,4λ,6λ-[1,3,5,2,4,6]三聚磷腈	2,2,4,4,6-Pentafluoro-6-phenoxy-2λ,4λ,6λ-[1,3,5,2,4,6]triazatriphosphinine	33027-68-8
2,3,4,5-四甲基-2-环戊烯酮	2,3,4,5-Tetramethyl-2-cyclopentenone	54458-61-6
二苯基二氟硅烷	DIPHENYLDIFLUOROSILANE	312-40-3
乙氧基(五氟)环三磷腈	Ethoxy(pentafluoro)cyclotriphosphazene	33027-66-6
全氟丁酸酐	Heptafluorobutyric anhydride	336-59-4
六氟环三磷腈	Hexafluorocyclotriphosphazene	15599-91-4
双氟磺酰亚胺锂盐 (LiFSI)	Lithium bis(fluorosulfonyl)imide	171611-11-3
油醇 (90N) 动物	Oleyl alcohol	143-28-2
苯氧基环磷腈	Phenoxy cyclophosphazene	1184-10-7

中文名称	英文名称	CAS
六氯三聚磷腈	PHOSPHONITRILIC CHLORIDE TRIMER	940-71-6
1,3-丙烯磺酸内酯 (PST)	Prop-1-ene-1,3-sultone	21806-61-1
四氟草酸磷酸锂	Tetrafluoroethylene oxalic acid phosphate lithium	521065-36-1
(OC-6-33')-Bis[6-(1H-pyrazol-1-yl-κN2)-2,2'-bipyridine-κN1,κN1']-Cobalt(2+)	(OC-6-33')-Bis[6-(1H-pyrazol-1-yl-κN2)-2,2'-bipyridine-κN1,κN1']-Cobalt(2+)	312322-08-0
(OC-6-33')-Bis[6-(1H-pyrazol-1-yl-κN2)-2,2'-bipyridine-κN1,κN1']-Cobalt(3+)	(OC-6-33')-Bis[6-(1H-pyrazol-1-yl-κN2)-2,2'-bipyridine-κN1,κN1']-Cobalt(3+)	1346417-60-4
2,3-萘酞菁铜		33273-09-5
2-[(7-{4-[N,N-Bis(4-methylphenyl)amino]phenyl}-2,1,3-benzothiadiazol-4-	2-[(7-{4-[N,N-Bis(4-methylphenyl)amino]phenyl}-2,1,3-benzothiadiazol-4-	1393343-58-2
2-[[7-(5-N,N-Ditolylaminothiophen-2-yl)-2,1,3-benzothiadiazol-4-yl]methylene]malononitrile	2-[[7-(5-N,N-Ditolylaminothiophen-2-yl)-2,1,3-benzothiadiazol-4-yl]methylene]malononitrile	1335150-09-8
N,N-二(9-菲基)-N,N-二苯基联苯胺		182507-83-1
Poly[(5,6-dihydro-5-octyl-4,6-dioxo-4H-thieno[3,4-c]pyrrole-1,3-diyl)[4,8-bis[(2-ethylhexyl)oxy]benzo[1,2-b:4,5-b']dithiophene-2,6-SMDPPO	Poly[(5,6-dihydro-5-octyl-4,6-dioxo-4H-thieno[3,4-c]pyrrole-1,3-diyl)[4,8-bis[(2-ethylhexyl)oxy]benzo[1,2-b:4,5-b']dithiophene-2,6-SMDPPO	1223479-75-1
聚(3,4-乙炔二氧噻吩)-聚(苯乙烯磺酸酯)		
聚(9,9-二正己基芴基-2,7-二基)		201807-75-2
聚[4,5-二氟-2,2-双(三氟甲基)-1,3-二氧杂环戊二烯-co-四氟乙炔]		37626-13-4
聚[9,9-双-(2-乙基己基)-9H-芴-2,7-二基]		188201-16-3
聚苯胺(翠绿亚胺盐)短链, 接枝聚合到木质素		313949-90-5
聚吡咯		30604-81-0
试剂DTS(PTTh2)2	DTS(PTTh2)2	1351374-94-1
试剂FBR	FBR	1644381-95-2
试剂FK 102 Co(II) PF6 salt	FK 102 Co(II) PF6 salt	1392221-69-0
试剂J61	J61	1887136-03-9
试剂LT-N875	LT-N875	939430-31-6
试剂PBDTTT-CF	PBDTTT-CF	1379592-65-0
试剂Poly(3,4-ethylenedioxythiophene), bis-poly(ethyleneglycol), lauryl terminated	Poly(3,4-ethylenedioxythiophene), bis-poly(ethyleneglycol), lauryl terminated	
试剂Polypyrrole-block-poly(caprolactone)	Polypyrrole-block-poly(caprolactone)	
试剂PTB7	PTB7	1266549-31-8
(5'-(3-(二苯基磷酰基)苯基)-[1,1':3',1"-三联苯]-3,3"-二基)双(二苯基氧化膦)	(5'-(3-(Diphenylphosphoryl)Phenyl)-[1,1':3',1"-Terphenyl]-3,3"-Diyl)Bis(Diphenylphosphine	1311378-95-6
3-氟-4-氧代吡啶-1-甲酸叔丁酯	[5,15-Bis(phenylethynyl)-10,20-bis[(triisopropylsilyl)ethynyl]porphyrinato]magnesium(II)	1397288-30-0
1,3,5-三(4-氟联苯-4-基)苯	1,3,5-Tris(4'-fluorobiphenyl-4-yl)benzene	372956-40-6
1,3,5-三[4-(9-咪唑基)苯基]苯	1,3,5-Tris[4-(9-Carbazolyl)Phenyl]Benzene	160780-82-5
1,3,5-三[4-(二苯基氨基)苯基]苯	1,3,5-Tris[4-(Diphenylamino)Phenyl]Benzene	147951-36-8
1,4,5,8-萘四甲酸酐	1,4,5,8-Naphthalenetetracarboxylic Dianhydride	81-30-1
1,4-二氨基丁烷二氢溴酸盐	1,4-Diaminobutane Dihydrobromide	18773-04-1
10,15-二氢-5,5,10,10,15,15-六甲基-5H-三苯并[a,f,k]三茚	10,15-Dihydro-5,5,10,10,15,15-hexamethyl-5H-tribenzo[a,f,k]trindene	597554-76-2
2,4-双[4-(二乙氨基)-2-羟基苯基]方酸	2-(4-(Diethylamino)-2-Hydroxyphenyl)-4-(4-(Diethylimino)-2-Hydroxycyclohexa-2,5-Dien-1-Ylidene)-3-Oxocyclobut-1-Enolate	68842-66-0
2-(4-叔丁苯基)-5-(4-联苯基)-1,3,4-恶二唑	2-(4-Tert-Butylphenyl)-5-(4-Biphenyl)-1,3,4-Oxadiazole	15082-28-7

中文名称	英文名称	CAS
9,9-二甲基-N,N'-双(3-甲基苯基)-N,N'-二苯基-9H-芴-2,7-二胺	2-(N-Methylaminocarbonyl)-1629-Pyridineboronic Acid Pincol Ester	677350-83-3
N,N'-双(萘-1-基)-N,N'-双(苯基)-2,7-二氨基-9,9-二甲基芴	2-(N-Methylaminocarbonyl)-1633-Pyridineboronic Acid Pincol Ester	1229226-27-0
N,N'-双(3-甲基苯基)-N,N'-双(苯基)-2,7-二氨基-9,9-二苯基芴	2-(N-Methylaminocarbonyl)-1637-Pyridineboronic Acid Pincol Ester	206886-03-5
4,4',4''-三(N-(萘-1-基)-N-苯基-氨基)三苯胺	2-(N-Methylaminocarbonyl)-1647-Pyridineboronic Acid Pincol Ester	185690-39-5
三(苯基吡唑)铱	2-(N-Methylaminocarbonyl)-1718-Pyridineboronic Acid Pincol Ester	562824-31-1
2,6-双(4-溴苯基)-4-(苊-1-基)-吡啶	2-(N-Methylaminocarbonyl)-1723-Pyridineboronic Acid Pincol Ester	
2,2',7,7'-四[N-萘基(苯基)-氨基]-9,9-螺双芴	2-(N-Methylaminocarbonyl)-1735-Pyridineboronic Acid Pincol Ester	404001-42-9
N,N'-双(菲-9-基)-N,N'-双(苯基)-联苯胺	2-(N-Methylaminocarbonyl)-1739-Pyridineboronic Acid Pincol Ester	934000-87-0
N,N,N',N'-四(2-萘基)联苯胺	2-(N-Methylaminocarbonyl)-1759-Pyridineboronic Acid Pincol Ester	186256-01-9
三(4-(喹啉-8-基)苯基)胺	2-(N-Methylaminocarbonyl)-1775-Pyridineboronic Acid Pincol Ester	1142945-07-0
N2,N7-二(萘-1-基)-9,9-二辛基-N2,N7-二苯基-9H-芴-2,7-二胺	2-(N-Methylaminocarbonyl)-1785-Pyridineboronic Acid Pincol Ester	870197-09-4
N,N'-双(3-甲基苯基)-N,N'-双(苯基)-9,9-二辛基芴	2-(N-Methylaminocarbonyl)-1789-Pyridineboronic Acid Pincol Ester	439942-97-9
N4,N4'-二(4-乙炔基苯基)-N4,N4'-二-1-萘基联苯-4,4'-二胺	2-(N-Methylaminocarbonyl)-1793-Pyridineboronic Acid Pincol Ester	1010396-31-2
N4,N4'-二(萘-1-基)-N4-(4-辛基苯基)-N4'-苯基联苯-4,4'-二胺	2-(N-Methylaminocarbonyl)-1797-Pyridineboronic Acid Pincol Ester	1431521-16-2
N4,N4'-双(4-(6-((3-乙基氧杂环丁-3-基)甲氧基)己基)苯基)-N4,N4'-二苯基联苯-4,4'-二胺	2-(N-Methylaminocarbonyl)-1801-Pyridineboronic Acid Pincol Ester	746634-00-4
N4,N4'-双(4-(6-((3-乙基氧杂环丁烷-3-基)甲氧基)己基)苯基)-N4,N4'-双(4-甲氧基苯基)-[1,1'-联	2-(N-Methylaminocarbonyl)-1805-Pyridineboronic Acid Pincol Ester	864130-79-0
4,4'-(二苯基亚甲基)双(N,N-二苯基苯胺)	2-(N-Methylaminocarbonyl)-1809-Pyridineboronic Acid Pincol Ester	1459723-98-8
N2,N7-二(萘-1-基)-N2,N7-二苯基-9,9-双(4-(4-乙炔基苄基)氧基)苯基)-9H-芴-2,7-二胺	2-(N-Methylaminocarbonyl)-1813-Pyridineboronic Acid Pincol Ester	1173170-48-3
N,N'-(4,4'-(环己烷-1,1-二基)双(4,1-亚苯基))双(N-(4-(6-(2-乙基氧杂环丁烷-2-基氧基)己基))苯基)-3,4,5-三氟苯胺	2-(N-Methylaminocarbonyl)-1817-Pyridineboronic Acid Pincol Ester	
聚[(9,9-二辛基芴基-2,7-二基)-alt-(N,N'-双{4-丁基苯基}-联苯胺-N,N'-{1,4-二苯基})]	2-(N-Methylaminocarbonyl)-1826-Pyridineboronic Acid Pincol Ester	223569-32-2
聚[(9,9-二己基芴基-2,7-二基)-alt-(N,N'-双-{p-丁基苯基}-1,4-二氨基苯基)]	2-(N-Methylaminocarbonyl)-1830-Pyridineboronic Acid Pincol Ester	870517-32-1
聚[(9,9-二辛基芴基-2,7-二基)-alt-(N,N'-双{对丁基苯基}-1,1'-联苯-4,4'-二胺)]	2-(N-Methylaminocarbonyl)-1834-Pyridineboronic Acid Pincol Ester	1115424-53-7
聚[2,7-[9,9-双(6'-(N,N-二乙氨基)己基)芴]-alt-2,7-[9,9-双(6-((3-乙基氧基-3-基)甲氧基)己基)]	2-(N-Methylaminocarbonyl)-1847-Pyridineboronic Acid Pincol Ester	1345045-27-3
2,2',7,7'-八(间甲苯胺)-9,9-螺二芴	2-(N-Methylaminocarbonyl)-1853-Pyridineboronic Acid Pincol Ester	302344-41-8
聚(9,9-双(6-三甲基碘化铵)-己基芴-2,7-二基-alt-(苯并[2,1,3]噻二唑-4,7-二基))	2-(N-Methylaminocarbonyl)-1857-Pyridineboronic Acid Pincol Ester	

中文名称	英文名称	CAS
4,4'-(二苯基硅烷二基)双(N,N-二苯基苯胺)	2-(N-Methylaminocarbonyl)-1861-Pyridineboronic Acid Pincol Ester	205327-13-5
4,4'-(9H-芴-9-亚烷基)双[N,N-双(4-甲基苯基)-苯甲胺]	2-(N-Methylaminocarbonyl)-1863-Pyridineboronic Acid Pincol Ester	159526-57-5
9,9-二甲基-10-(9-苯基-9H-吡啶-3-基)-9,10-二氢吡啶	2-(N-Methylaminocarbonyl)-1867-Pyridineboronic Acid Pincol Ester	1705584-08-2
2,2'-二甲基-N4,N4,N4',N4'-四间甲基联苯-4,4'-二胺	2-(N-Methylaminocarbonyl)-1871-Pyridineboronic Acid Pincol Ester	80730-98-9
聚[(9,9-二辛基芴基-2,7-二基)-co-(4,4'-(N-(2,4-二甲基苯基)二苯胺)]	2-(N-Methylaminocarbonyl)-1879-Pyridineboronic Acid Pincol Ester	
聚[(9,9-二辛基芴基-2,7-二基)-co-(4,4'-(N-(二苯胺)-(2-氰基异丙基苯基)]	2-(N-Methylaminocarbonyl)-1883-Pyridineboronic Acid Pincol Ester	
3,5-二-9H-吡啶-9-基-N,N-双[4-[[6-[(3-乙基-3-氧代坦基)甲氧基]己基]氧基]苯基]苯胺	2-(N-Methylaminocarbonyl)-1887-Pyridineboronic Acid Pincol Ester	2351818-31-8
N4,N4'-(联苯-4,4'-二基)双(N4'-(萘-1-基)-N4,N4'-二苯基联苯-4,4'-二胺)	2-(N-Methylaminocarbonyl)-1891-Pyridineboronic Acid Pincol Ester	292827-46-4
N2,N2'-(9,9-二甲基-9H-芴-2,7-二基)双(9,9-二甲基-N2,N7,N7'-三苯基-9H-芴-2,7-二胺)	2-(N-Methylaminocarbonyl)-1895-Pyridineboronic Acid Pincol Ester	354987-71-6
N1,N1'-(联苯基-4,4'-二基)双(N1-(萘-2-基)-N4,N4'-二苯基苯-1,4-二胺)	2-(N-Methylaminocarbonyl)-1903-Pyridineboronic Acid Pincol Ester	1207378-72-0
双吡嗪[2,3-f:2',3'-h]喹啉啉-2,3,6,7,10,11-六碳腈	2-(N-Methylaminocarbonyl)-1907-Pyridineboronic Acid Pincol Ester	215611-93-1
2,2'-双[N,N-双(4-甲氧基苯基)氨基]-9,9-螺双芴	2-(N-Methylaminocarbonyl)-1919-Pyridineboronic Acid Pincol Ester	1174006-40-6
N,N'-二苯基-N,N'-二-[4-(N,N-二对甲苯基氨基)苯基]联苯胺	2-(N-Methylaminocarbonyl)-1926-Pyridineboronic Acid Pincol Ester	1394130-64-3
N,N'-二苯基-N,N'-二-[4-(N,N-二苯基氨基)苯基]联苯胺	2-(N-Methylaminocarbonyl)-1927-Pyridineboronic Acid Pincol Ester	936355-01-0
2,3,8,9,14,15-六氯二喹啉并[2,3-a:2',3'-c]吩嗪	2-(N-Methylaminocarbonyl)-1947-Pyridineboronic Acid Pincol Ester	389121-44-2
2,3,8,9,14,15-六氟-5,6,11,12,17,18-六氮杂联三萘	2-(N-Methylaminocarbonyl)-1951-Pyridineboronic Acid Pincol Ester	872140-95-9
N2,N2'-(9,9-二辛基-9H-芴-2,7-二基)双(9,9-二甲基-N2,N7,N7'-三苯基-9H-芴-2,7-二胺)	2-(N-Methylaminocarbonyl)-1955-Pyridineboronic Acid Pincol Ester	
2-(2-甲氧基苯基)-1,3-二甲基-1H-苯并咪唑-3-碘鎓	2-(N-Methylaminocarbonyl)-1959-Pyridineboronic Acid Pincol Ester	
4,7-二苯基-2,9-双(4-(1-苯基-1H-苯并[d]咪唑-2-基)苯基)-1,10-菲咯啉	2-(N-Methylaminocarbonyl)-3007-Pyridineboronic Acid Pincol Ester	1447848-17-0
2,4,6-三(3-(咪啉-5-基)苯基)-1,3,5-三嗪	2-(N-Methylaminocarbonyl)-3011-Pyridineboronic Acid Pincol Ester	1874199-82-2
聚{[(9,9-双[3-(三甲基铵)丙基-2,7-芴]-alt-2,7-(9,9-二辛基芴)]二碘化物}	2-(N-Methylaminocarbonyl)-3015-Pyridineboronic Acid Pincol Ester	1687752-60-8
聚[(9,9-双(3'-(N,N-三甲基氨基)丙基)-2,7-芴)-alt-1,4-苯]二碘化物	2-(N-Methylaminocarbonyl)-3019-Pyridineboronic Acid Pincol Ester	1026670-72-3
4,6-双(3,5-二(2-吡啶)基苯基)-2-甲基咪啉	2-(N-Methylaminocarbonyl)-3023-Pyridineboronic Acid Pincol Ester	1266181-51-4
1,3-双[2-(2,2'-联吡啶-6-基)-1,3,4-恶二唑-5-基]苯	2-(N-Methylaminocarbonyl)-3027-Pyridineboronic Acid Pincol Ester	866117-19-3
6,6'-双[5-(联苯-4-基)-1,3,4-恶二唑-2-基]-2,2'-联吡啶基	2-(N-Methylaminocarbonyl)-3031-Pyridineboronic Acid Pincol Ester	1219827-28-7



中文名称	英文名称	CAS
6-甲基-2-(4-(9-(4-(6-甲基苯并[d]噻唑-2-基)苯基)噻-10-基)苯基)苯并[d]噻唑	2-(N-Methylaminocarbonyl)-3035-Pyridineboronic Acid Pincol Ester	850018-19-8
2,7-双[2-(2,2'-联吡啶-6-基)-1,3,4-恶二唑-5-基]-9,9-二甲基芴	2-(N-Methylaminocarbonyl)-3044-Pyridineboronic Acid Pincol Ester	1174006-45-1
2-(萘-2-基)-4,7-二苯基-1,10-菲咯啉	2-(N-Methylaminocarbonyl)-3050-Pyridineboronic Acid Pincol Ester	923972-84-3
苯基二-1-苊基氧磷	2-(N-Methylaminocarbonyl)-3054-Pyridineboronic Acid Pincol Ester	721969-93-3
3,3',5,5'-四[(间吡啶基)-苯-3-基]联苯	2-(N-Methylaminocarbonyl)-3058-Pyridineboronic Acid Pincol Ester	1009033-94-6
二苯基双(4-(吡啶-3-基)苯基)硅烷	2-(N-Methylaminocarbonyl)-3068-Pyridineboronic Acid Pincol Ester	1152162-74-7
1,3,5-三(4-(吡啶-4-基)喹啉-2-基)苯	2-(N-Methylaminocarbonyl)-3081-Pyridineboronic Acid Pincol Ester	1350742-68-5
2-(2',2''-联吡啶-6'-基)酚酸锂	2-(N-Methylaminocarbonyl)-3092-Pyridineboronic Acid Pincol Ester	1049805-81-3
4,6-双(3,5-二(吡啶-4-基)苯基)-2-苯基嘧啶	2-(N-Methylaminocarbonyl)-3096-Pyridineboronic Acid Pincol Ester	1097652-83-9
1,3,5-三(6-(3-(吡啶-3-基)苯基)吡啶-2-基)苯	2-(N-Methylaminocarbonyl)-3097-Pyridineboronic Acid Pincol Ester	1492917-78-8
4,6-双(3,5-二(吡啶-3-基)苯基)-2-(吡啶-3-基)嘧啶	2-(N-Methylaminocarbonyl)-3099-Pyridineboronic Acid Pincol Ester	1382639-67-9
4,6-双(3,5-二(吡啶-4-基)苯基)-2-(3-(吡啶-3-基)苯基)嘧啶	2-(N-Methylaminocarbonyl)-3100-Pyridineboronic Acid Pincol Ester	1382639-70-4
聚[9,9-双[6'-(N,N,N-三甲基铵)己基]芴-alt-co-1,4-亚苯基]溴化物	2-(N-Methylaminocarbonyl)-3108-Pyridineboronic Acid Pincol Ester	
8-羟基喹啉钠盐	2-(N-Methylaminocarbonyl)-3112-Pyridineboronic Acid Pincol Ester	2872-54-0
聚[(9,9-双(2-乙基己基)芴基-2,7-二基)-alt-(苯并[2,1,3]噻二唑-4,7-二基)]	2-(N-Methylaminocarbonyl)-3116-Pyridineboronic Acid Pincol Ester	1379051-13-4
聚(9,9-二正辛基芴基-2,7-亚苯基乙撑)	2-(N-Methylaminocarbonyl)-3124-Pyridineboronic Acid Pincol Ester	360068-36-6
聚[N-9'-庚烷基-2,7-吡啶-alt-3,6-双(噻吩-5-基)-2,5-二辛基-2,5-二氢吡咯[3,4]吡咯-1,4-二酮]	2-(N-Methylaminocarbonyl)-3152-Pyridineboronic Acid Pincol Ester	1147374-60-4
(E)-4-((5-羧基-3,3-二甲基-1-辛基-3H-吡啶-2-基)亚甲基)-2-((E)-(1-乙基-3,3-二甲基吡啶-2-亚基)甲基)-3-氧代环己烷-1-烯酸盐	2-(N-Methylaminocarbonyl)-3164-Pyridineboronic Acid Pincol Ester	951248-19-4
(E)-4-((5-羧基-3,3-二甲基-1-辛基-3H-吡啶鎓-2-基)亚甲基)-2-((E)-(3-乙基-1,1-二甲基-1H-苯并[e]吡啶-2(3H)-亚基)甲基)-3-氧代环丁-1-烯醇	2-(N-Methylaminocarbonyl)-3168-Pyridineboronic Acid Pincol Ester	1240041-84-2
3,6-双(5-(苯并咪唑-2-基)噻吩-2-基)-2,5-双(2-乙基己基)吡咯[3,4-c]吡咯-1,4(2H,5H)-二酮	2-(N-Methylaminocarbonyl)-3176-Pyridineboronic Acid Pincol Ester	1241836-03-2
2-((Z)-2-((E)-2-(1,1-二甲基-5,6-二氢-1H-吡咯并[3,2,1-IJ]喹啉-2(4H)-亚烷基)亚乙基)-3-氧代-2,3-二氢-1H-茛-1-基)丙二腈	2-(N-Methylaminocarbonyl)-3200-Pyridineboronic Acid Pincol Ester	1207961-82-7
2-(6-(4-(二苯氨基)苯基)-4,4'-二苯基二硫代[3,2-b:2',3'-d]硅氧基-2-基)-亚甲基丙二腈	2-(N-Methylaminocarbonyl)-3217-Pyridineboronic Acid Pincol Ester	1321562-79-1
2,2'-(5,5'-(噻吩并[3,2-b]噻吩-2,5-二基)双(噻吩-5,2-二基))双(甲烷-1-基-1-亚烷基)二马腈	2-(N-Methylaminocarbonyl)-3237-Pyridineboronic Acid Pincol Ester	1432613-96-1
聚[2,6-(4,4'-双(乙基己基)二硫杂环[3,2-b:2',3'-d]硅油)-alt-(1,3-(5-辛基-4H-噻吩[3,4-c]吡咯-	2-(N-Methylaminocarbonyl)-3241-Pyridineboronic Acid Pincol Ester	1279109-93-1

中文名称	英文名称	CAS
2-氰基-3-(7-(4-(二苯基氨基)苯基)-苯并[C][1,2,5]噻二唑-4-基)丙烯酸	2-(N-Methylaminocarbonyl)-3245-Pyridineboronic Acid Pincol Ester	1415134-59-6
2,4-双(4-(二异丁基氨基)-2,6-二羟基苯基)环丁烷-1,3-双(酚盐)	2-(N-Methylaminocarbonyl)-3253-Pyridineboronic Acid Pincol Ester	432493-75-9
N4,N4'-双(二苯并[b,d]噻吩-4-基)-N4,N4'-二苯基联苯-4,4'-二胺	2-(N-Methylaminocarbonyl)-3277-Pyridineboronic Acid Pincol Ester	1092536-44-1
2-(3,6-双(二-对甲苯胺基)-9H-芴-9-亚烷基)丙二腈	2-(N-Methylaminocarbonyl)-3281-Pyridineboronic Acid Pincol Ester	1527494-19-4
聚{2,2'-[(2,5-二(2-己基癸基)-3,6-二氧基-2,3,5,6-四氢吡咯[3,4-c]吡咯-1,4-二酰基)二硫代酚]-5,5'-二酰基-alt-噻吩-2,5-二酰基}	2-(N-Methylaminocarbonyl)-3338-Pyridineboronic Acid Pincol Ester	1198291-01-8
聚[3,6-双(4'-十二烷基[2,2'-联噻吩]-5-基)-2,5-双(2-乙基己基)-2,5-二氢吡咯并[3,4-c]-吡咯-1,4-二	2-(N-Methylaminocarbonyl)-3350-Pyridineboronic Acid Pincol Ester	1002094-22-5
聚{[2,7-(9,9-二(3,7-二甲基辛基芴))-alt-[5,5-(4,7-二-2'-噻吩基-2,1,3-苯并噻唑)]	2-(N-Methylaminocarbonyl)-3354-Pyridineboronic Acid Pincol Ester	1204005-82-2
聚{2,2'-[(2,5-双(2-己基癸基)-3,6-二氧基-2,3,5,6-四氢吡咯[3,4-c]吡咯-1,4-二基)]二硫苯-5,5'-二酰基-alt-苯-1,4-二基}	2-(N-Methylaminocarbonyl)-3358-Pyridineboronic Acid Pincol Ester	1255939-39-9
聚{2,2'-[(2,5-二(2-辛基十二烷基)-3,6-二氧基-2,3,5,6-四氢吡咯[3,4-c]吡咯-1,4-二基)]二硫苯-5,5'-二酰基噻吩[3,2-b]噻吩-2,5-二基}	2-(N-Methylaminocarbonyl)-3366-Pyridineboronic Acid Pincol Ester	1260685-65-1
聚并二噻吩吡咯并吡咯二酮-噻吩	2-(N-Methylaminocarbonyl)-3370-Pyridineboronic Acid Pincol Ester	1270977-98-4
2,6-双[2,5-双(3-辛基罗丹宁)-(3,3-二辛基-2,2:5,2'-对噻吩)]-4,8-双((5-乙基己基)噻吩-2-基)苯并[1,2-B:4,5-B]二噻吩	2-(N-Methylaminocarbonyl)-3378-Pyridineboronic Acid Pincol Ester	1651168-29-4
3',4'-二丁基-5,5''-双(二氰基乙烯基)-2,2':5',2''-对噻吩	2-(N-Methylaminocarbonyl)-3386-Pyridineboronic Acid Pincol Ester	908588-68-1
三[4-叔丁基-2-(1H-吡唑-1-基)吡啶]钴三(1,1,1-三氟-N-[(三氟甲基)磺酰基]甲烷磺酰胺盐)	2-(N-Methylaminocarbonyl)-3391-Pyridineboronic Acid Pincol Ester	1447938-61-5
N2,N2',N7,N7'-四(3-甲氧基苯基)-N2,N2',N7,N7'-四(4-甲氧基苯基)-9,9'-螺双[芴]-2,2',7,7'-四胺	2-(N-Methylaminocarbonyl)-3425-Pyridineboronic Acid Pincol Ester	1573202-44-4
N2,N2',N7,N7'-四(2-甲氧基苯基)-N2,N2',N7,N7'-四(4-甲氧基苯基)-9,9'-螺双[芴]-2,2',7,7'-四胺	2-(N-Methylaminocarbonyl)-3429-Pyridineboronic Acid Pincol Ester	1628961-22-7
2,3,9,1,16,17-六氯化硼次酞菁硼	2-(N-Methylaminocarbonyl)-3433-Pyridineboronic Acid Pincol Ester	309963-65-3
2,2',7,7'-四(4,4'-二甲氧基-3-甲基二苯基氨基)-9,9'-螺二芴	2-(N-Methylaminocarbonyl)-3453-Pyridineboronic Acid Pincol Ester	1573202-31-9
锌(II)-1,8,15,22-四氟-29H,31H-酞菁,异构体混合物	2-(N-Methylaminocarbonyl)-3457-Pyridineboronic Acid Pincol Ester	1120355-28-3
锌(II)-2,3,9,1,16,17,23,24-八氟-29H,31H-酞菁	2-(N-Methylaminocarbonyl)-3461-Pyridineboronic Acid Pincol Ester	676519-80-5
铁(II)1,2,3,4,8,9,1,11,15,16,17,18,22,23,24,25-十六氟-29H,31H-酞菁	2-(N-Methylaminocarbonyl)-3477-Pyridineboronic Acid Pincol Ester	23844-93-1
聚[[2,5-二(2-己基癸基)-2,3,5,6-四氧-3,6-二氧吡咯[3,4-c]吡咯-1,4-二酰基]-alt-[3',3'-二甲基-2,2':5',2'-三噻吩]-5,5'-二酰基]	2-(N-Methylaminocarbonyl)-3481-Pyridineboronic Acid Pincol Ester	1429119-71-0
2,6-双[(5'-辛基[2,2'-联噻吩]){4,8-双(5-辛基-2-噻吩基)苯并[1,2-b:4,5-b']二噻吩]	2-(N-Methylaminocarbonyl)-3485-Pyridineboronic Acid Pincol Ester	

中文名称	英文名称	CAS
[2,3,9,10,16,17-六氯-7,12:14,19-二亚氨基-21,5-硝基-5H-三苯并[c,h,m] [1,6,11]三氮杂环戊四烯酸酯(2)-N22,N23,N24]苯氧基-硼	2-(N-Methylaminocarbonyl)-3493-Pyridineboronic Acid Pincol Ester	1309390-01-9
三(2,2'-联吡啶)双(三氟甲磺酰亚胺)钴(II)	2-(N-Methylaminocarbonyl)-3509-Pyridineboronic Acid Pincol Ester	1353745-90-0
三(2,2'-联吡啶)三(三氟甲磺酰亚胺)钴(III)	2-(N-Methylaminocarbonyl)-3513-Pyridineboronic Acid Pincol Ester	1353745-91-1
6'',5'',10'':4',5',6'二异喹啉并[2,1-A:2',1'-A']二萘嵌间二氮杂苯-12,25-二酮	2-(N-Methylaminocarbonyl)-3525-Pyridineboronic Acid Pincol Ester	6859-32-1
三(5-硝基-1,10-菲咯啉)双(六氟磷酸盐)钴(II)	2-(N-Methylaminocarbonyl)-3537-Pyridineboronic Acid Pincol Ester	
三(5-硝基-1,10-菲咯啉)三(六氟磷酸盐)三(六氟磷酸盐)	2-(N-Methylaminocarbonyl)-3541-Pyridineboronic Acid Pincol Ester	
4-甲氧基苯基碘化铵	2-(N-Methylaminocarbonyl)-3545-Pyridineboronic Acid Pincol Ester	
二异丙基碘化铵	2-(N-Methylaminocarbonyl)-3546-Pyridineboronic Acid Pincol Ester	
己胺氢碘酸盐	2-(N-Methylaminocarbonyl)-3547-Pyridineboronic Acid Pincol Ester	54285-91-5
3,4,9,10-Per四甲酸二酰亚胺	2-(N-Methylaminocarbonyl)-3548-Pyridineboronic Acid Pincol Ester	891911-18-5
六氟磷酸甲脒	2-(N-Methylaminocarbonyl)-3571-Pyridineboronic Acid Pincol Ester	
四氟硼酸甲脒	2-(N-Methylaminocarbonyl)-3572-Pyridineboronic Acid Pincol Ester	
N,N'-双(N-庚基)-3,4,9,10-二甲酰亚胺	2-(N-Methylaminocarbonyl)-3589-Pyridineboronic Acid Pincol Ester	95689-91-1
聚[萘酚[1,2-c:5,6-c']双[1,2,5]噻二唑-5,10-二烷基[3,3''-双(2-癸基十四烷基)[2,2':5',2':5',2'-四噻吩]-	2-(N-Methylaminocarbonyl)-3625-Pyridineboronic Acid Pincol Ester	1357999-94-0
聚[4,8-双{(N,N-二甲基)-N-乙基铵}丙氧基]苯并[1,2-b:4,5-b]二硫代苯-2,6-二酰基-alt-3-氟-2-[(2-乙基己基)羰基]噻吩[3,4-b]噻吩-4,6-二酰基]	2-(N-Methylaminocarbonyl)-3633-Pyridineboronic Acid Pincol Ester	2092936-05-3
2,5,8-三(4-(双(4-甲氧基苯基)氨基)苯基)苯并三噻吩	2-(N-Methylaminocarbonyl)-3637-Pyridineboronic Acid Pincol Ester	1908506-93-3
聚[[5,6-二(十二烷氧基)-2,1,3-苯并噻唑-4,7-二基]-2,5-噻吩二基[4,8-二[4-(2-乙基己基)-3,5-二氟苯基]苯并[1,2-b:4,5-b]二噻吩-2,6-二基]-2,5-	2-(N-Methylaminocarbonyl)-3641-Pyridineboronic Acid Pincol Ester	1840869-86-4
氧代酞菁锡(IV)	2-(N-Methylaminocarbonyl)-3661-Pyridineboronic Acid Pincol Ester	57156-41-9
聚[2,2''-双[(2-丁基氧基)羰基][2,2':5',2''-季硫酚]-5,5''-二酰基]	2-(N-Methylaminocarbonyl)-3671-Pyridineboronic Acid Pincol Ester	1609536-17-5
聚[[2,7-二(2-乙基己基)-1,2,3,6,7,8-六氢-1,3,6,8-四氧代苯并[Imn][3,8]邻菲罗啉-4,9-二酰基]-2,5-噻吩二酰基[9,9-二[3'((N,N-二甲基)-N-乙基铵)-丙基]-9H-氟苯-2,7-二酰基]-2,5-噻吩二酰基]	2-(N-Methylaminocarbonyl)-3675-Pyridineboronic Acid Pincol Ester	2169941-79-9
4,4,4,4-硅烷四基四(N,N-双(4-甲氧基苯基)苯胺)	2-(N-Methylaminocarbonyl)-3679-Pyridineboronic Acid Pincol Ester	
新戊基碘化铵	2-(N-Methylaminocarbonyl)-3687-Pyridineboronic Acid Pincol Ester	
4-氟苄基溴化铵	2-(N-Methylaminocarbonyl)-3690-Pyridineboronic Acid Pincol Ester	

中文名称	英文名称	CAS
聚(9,9-二辛基-2,7-二基)	2-(N-Methylaminocarbonyl)-3692-Pyridineboronic Acid Pincol Ester	123864-00-6
二甲基苯基封端的聚(9,9-二辛基氟苯基-2,7-二酰基)	2-(N-Methylaminocarbonyl)-3696-Pyridineboronic Acid Pincol Ester	874816-14-5
4-三氟甲基苯基溴化铵	2-(N-Methylaminocarbonyl)-3697-Pyridineboronic Acid Pincol Ester	
1,4-二氮杂双环[2,2,2]辛烷-1,4-二四氟硼酸盐	2-(N-Methylaminocarbonyl)-3708-Pyridineboronic Acid Pincol Ester	
四氟硼酸胍	2-(N-Methylaminocarbonyl)-3713-Pyridineboronic Acid Pincol Ester	36595-00-3
聚[[2,2'-双榄香烯]-5,5'-二基(9,9-二辛基-9H-苈-2,7-二基)]	2-(N-Methylaminocarbonyl)-3736-Pyridineboronic Acid Pincol Ester	898233-56-2
聚[(噻吩)-alt-(6,7-二氟-2-(2-己基癸氧基)喹啉)]	2-(N-Methylaminocarbonyl)-3780-Pyridineboronic Acid Pincol Ester	2270233-86-6
N2,N2',N7,N7'-四(2,4-二甲氧基苯基)-N2,N2',N7,N7'-四苯基-9,9'-螺环[苈]-2,2',7,7'-四	2-(N-Methylaminocarbonyl)-3784-Pyridineboronic Acid Pincol Ester	2102348-92-3
N,N-二苯基-[1,1'-联苯]-4-胺-4'-[4,5-二苯基-1-[4-(10-苯基-9-4'-[4,5-二苯基-1-[4-(10-苯基-9-萸基)苯基]-1H-咪唑-2-基]	2-(N-Methylaminocarbonyl)-3788-Pyridineboronic Acid Pincol Ester	2252254-01-4
聚[[1,2,3,6,7,8-六氢-2,7-双(2-辛基十二烷基)-1,3,6,8-四氧代苯并[Imn][3,8]邻菲罗啉-4,9-二基]-[3,3'-二氟[2,2'-双噻吩]-5,5'-二基]]	2-(N-Methylaminocarbonyl)-3793-Pyridineboronic Acid Pincol Ester	1773524-46-1
聚{[N,N'-双(2-己基癸基)萘-1,4,5,8-双(二甲酰亚胺)-2,6-二酰基]-alt-2,5-噻吩}	2-(N-Methylaminocarbonyl)-3805-Pyridineboronic Acid Pincol Ester	1675973-74-6
聚[[1,2,3,6,7,8-六氢-2,7-二(2-辛基十二烷基)-[1,3,6,8-四氧代苯并[Imn][3,8]菲罗啉-4,9-二基]-	2-(N-Methylaminocarbonyl)-3809-Pyridineboronic Acid Pincol Ester	1236111-76-4
聚[(5,6-二氢-5-辛基-4,6-二氧基-4H-噻吩[3,4-c]吡咯-1,3-二酰基)[4,8-双[5-(2-乙基己基,2-噻吩基)苯并[1,2-b:4,5-b']二硫代-2,6-二酰基]	2-(N-Methylaminocarbonyl)-3813-Pyridineboronic Acid Pincol Ester	1426534-44-2
聚[[4,8-二[4-氯-5-(2-乙基己基和2-噻吩基)苯并[1,2-b:4,5-b']二硫代苯-2,6-二基]-2,5-噻吩二基[5,7-二(2-乙基己基,4,8-二氧基-4H,8H-苯并[1,2-c:4,5-c']二硫代苯-1,3-二基)-2,5-噻吩二基]	2-(N-Methylaminocarbonyl)-3818-Pyridineboronic Acid Pincol Ester	2239295-71-5
聚[[2,5-二(2-癸基十四烷基)-2,3,5,6-四氢-3,6-二氧吡咯[3,4-c]吡咯-1,4-二酰基](3",4'-二氟[2,2':5",2":5",2'-四氢苯基]-5'-二酰基)]	2-(N-Methylaminocarbonyl)-3827-Pyridineboronic Acid Pincol Ester	1926166-46-2
聚[噻唑并[5,4-d]噻唑-2,5-二基[4-(2-乙基己基)-2,5-噻吩二酰基][4,8-双[5-(2-丁基辛基)-2-噻吩基]-苯并[1,2-b:4,5-b']二噻吩-2,6-二基][3-(2-乙基己基)-2,5-噻吩二酰基]]	2-(N-Methylaminocarbonyl)-3831-Pyridineboronic Acid Pincol Ester	2029196-34-5
(5Z,5'Z)-5,5'-((7,7'-(6,6,12,12-四辛基-6,12-二氢茛并[1,2-b]苈-2,8-二酰基)双(苯并[c][1,2,5]噻二唑-7,4-二酰基)双(甲酰亚砷)双(3-乙基-2-硫氧噻唑	2-(N-Methylaminocarbonyl)-3839-Pyridineboronic Acid Pincol Ester	2088628-64-0
聚[[5,6-二氟-2-(2-己基癸基)-2H-苯并三唑-4,7-二基]-2,5-噻吩二基[4,8-双[5-(2-乙基己基)-2-噻吩基]苯并[1,2-b:4,5-b']二噻吩-2,6-二基]-2,5-噻吩	2-(N-Methylaminocarbonyl)-3843-Pyridineboronic Acid Pincol Ester	1887136-01-7
硼亚2,3-萘酞菁氯化物	2-(N-Methylaminocarbonyl)-3847-Pyridineboronic Acid Pincol Ester	142710-56-3
3,9-双(2-亚甲基-((3-(1,1-二氟基亚甲基)-6/7-甲基)-茛满酮))-5,5,11,11-四(4-己基苯基)-二噻吩[2,3-d:2',3'-d']-s-茛三烯[1,2-b:5,6-b']二噻吩	2-(N-Methylaminocarbonyl)-3855-Pyridineboronic Acid Pincol Ester	

中文名称	英文名称	CAS
2,2'-((2Z,2'Z)-((4,4,4',4'-四己基-4H,4'H-[2,2'-bi(环戊二烯[1,2-b:5,4-b']二噻吩))-6,6'-二基)双(亚甲基))双(5,6-二氟-3-氧代-2,3-二氢-1H-茛-2,1-二亚	2-(N-Methylaminocarbonyl)-3859-Pyridineboronic Acid Pincol Ester	
2,9-双(1,2-二氢茛-5-基)-5,12-二辛基-5,12-二氢喹啉[2,3-b]吡啶-7,14-二酮	2-(N-Methylaminocarbonyl)-3872-Pyridineboronic Acid Pincol Ester	
2,2'-((2Z,2'Z)-((12,13-二(2-乙基己基)-3,9-二十一烷基-12,13-二氢-[1,2,5]噻二唑[3,4-e]噻吩-[2,3':4',5']-噻吩[2',3':4,5]吡咯[3,2-g]噻吩-[2',3':4,5]噻吩[3',-b]-吡啶-2,10-二酰基)双(甲酰亚砷)-双(5,6-二氯-3-3-3-二氢-1,1,2'-2'-1-二氢吡	2-(N-Methylaminocarbonyl)-3880-Pyridineboronic Acid Pincol Ester	2447642-40-0
3,9-双(2-亚甲基-(3-(1,1-二氧基亚甲基)-6/7-氟)-茛满酮))-5,5,11,11-四(5-己基噻吩基)-二噻吩[2,3-d:2',3'-d']-s-茛二烯[1,2-b:5,6-b']二噻吩	2-(N-Methylaminocarbonyl)-3884-Pyridineboronic Acid Pincol Ester	
3,9-双(2-亚甲基-(3-(1,1-二氧基亚甲基)-茛满酮))-5,5,11,11-四(3-己基苯基)-二噻吩并[2,3-d:2',3'-d']-s-茛三烯[1,2-b:5,6-b']-二噻吩	2-(N-Methylaminocarbonyl)-3888-Pyridineboronic Acid Pincol Ester	
试剂ITIC-Cl	2-(N-Methylaminocarbonyl)-3892-Pyridineboronic Acid Pincol Ester	2253663-81-7
3,9-双(2-亚甲基-(3-(1,1-二氧基亚甲基)-6/7-氟)-茛满酮))-5,5,11,11-四(4-己基苯基)-二噻吩[2,3-d:2',3'-d']-s-茛三烯[1,2-b:5,6-b']二噻吩	2-(N-Methylaminocarbonyl)-3896-Pyridineboronic Acid Pincol Ester	
3,9-双(2-亚甲基-(3-(1,1-二氧基亚甲基)-6/7-氟)-茛满酮))-5,5,11,11-四(4-己基苯基)-二噻吩[2,3-d:2',3'-d']-s-茛三烯[1,2-b:5,6-b']二噻吩	2-(N-Methylaminocarbonyl)-3900-Pyridineboronic Acid Pincol Ester	
2,2'-((2Z,2'Z)-((12,13-双(2-乙基己基)-3,9-二十一烷基-12,13-二氢-2-(2-乙基己基)-[1,2,3]三氮唑[3,4-e]噻吩并[2',3':4',5']噻吩并[2',3':4,5]吡咯并[3,2-g]噻吩并[2',3':4,5]噻吩并-[3,2-b]吡啶-2,10-二基)双(甲酰亚甲基)双-(5,6-二氟-3-氧代-2,3-二	2-(N-Methylaminocarbonyl)-3912-Pyridineboronic Acid Pincol Ester	
[聚[2,1,3-苯并噻二唑-4,7-二基[4,4-双(2-乙基己基)-4H-硅杂环戊二烯并[3,2-B4,5-B']二噻吩-2,6-	2-(N-Methylaminocarbonyl)-3916-Pyridineboronic Acid Pincol Ester	1089687-02-4
聚[[4,4-双(2-乙基己基)-4H-环戊并[2,1-B:3,4-B]二噻吩-2,6-二基]-2,5-噻吩二基-2,1,3-苯并噻二唑-4,7-二基-2,5-噻吩二基]	2-(N-Methylaminocarbonyl)-3924-Pyridineboronic Acid Pincol Ester	1035681-83-4
2,2'-((2Z,2'Z)-((12,13-双(2-乙基己基)-3,9-二十一烷基-12,13-二氢-2-(2-乙基己基)-[1,2,3]三氮唑[3,4-e]噻吩并[2',3':4',5']噻吩并[2',3':4,5]吡咯并[3,2-g]噻吩并[2',3':4,5]噻吩并-[3,2-b]吡啶-2,10-二基)-双(甲酰亚基)双-(5,6-二氟-3-氧代-2,3-二	2-(N-Methylaminocarbonyl)-3932-Pyridineboronic Acid Pincol Ester	
9-乙基-N2,N7-双(9-(2-乙基己基)-9H-咪唑-3-基)-N2,N7-双(4-甲氧基苯基)-9H-咪唑-2,7-二胺盐酸	2-(N-Methylaminocarbonyl)-3940-Pyridineboronic Acid Pincol Ester	
2,2-[[12,13-双(2-丁基辛基)-12,13-二氢-3,9-二壬基双噻吩并[2,3:4,5]噻吩并[2,3:4,5]吡咯[3,2-E:2,3-G][2,1,3]苯并噻二唑-2,10-二基]双[甲基亚甲基(5,6-氯-3-氧-1H-茛-2,1(3H)-二亚甲基)]双	2-(N-Methylaminocarbonyl)-3948-Pyridineboronic Acid Pincol Ester	2598965-39-8
N2,N2,N7,N7-四(3-氟-4-甲氧基苯基)-N2,N2,N7,N7-四(4-甲氧基苯基)-9,9'-茛-2,2',7,7'-2,5-二己基-3,6-双(5"-正己基-[2,2',5',2"]对噻吩-5-基)-吡咯并[3,4-C]吡咯-1,4-二酮	2-(N-Methylaminocarbonyl)-3953-Pyridineboronic Acid Pincol Ester	
	2-(N-Methylaminocarbonyl)-4099-Pyridineboronic Acid Pincol Ester	1077629-69-6

中文名称	英文名称	CAS
2,2'-((2Z, 2'Z)-((4,4,9,9)-四己基-4,9-二氢-s-茛并[1,2-b:5,6-b']二噻吩-2,7-二基)双(甲酰亚胺)双(3-氧代-2,3-二氢-1小时-茛-2,1-二亚基))二丙二腈	2,2'-((2Z,2'Z)-((4,4,9,9-Tetrahexyl-4,9-Dihydro-s-Indaceno[1,2-b:5,6-b']Dithiophene-2,7-Diyl)Bis(Methanylylidene))Bis(3-Oxo-2,3-Dihydro-1H-Indene-2,1-	1883441-92-6
试剂2,2'-[[6,6,12,12-tetrakis(4-hexylphenyl)-s-indacenodithieno[3,2-b]thiophene]methylidyne(3-oxo-1H-indene-2,1(3H)-	2,2'-[[6,6,12,12-tetrakis(4-hexylphenyl)-s-indacenodithieno[3,2-b]thiophene]methylidyne(3-oxo-1H-indene-2,1(3H)-	1664293-06-4
2,4,6-三([1,1'-联苯]-4-基)-1,3,5-三嗪	2,4,6-Tri([1,1'-Biphenyl]-4-Yl)-1,3,5-Triazine	31274-51-8
2,4,6-三(9H-咔唑-9-基)-1,3,5-三嗪	2,4,6-Tri(9H-Carbazol-9-Yl)-1,3,5-Triazine	134984-37-5
2,4,6-三苯基三嗪	2,4,6-Triphenyl-1,3,5-Triazine	493-77-6
2,4-双[8-羟基-1,1,7,7-四甲基久洛尼定-9-基]方酸	2,4-Bis[8-Hydroxy-1,1,7,7-Tetramethyljulolidin-9-Yl]Squaraine	358727-55-6
2,5-双(2-乙基己基)-3,6-双(5-己基[2,2:5,2-三联噻吩]-5-基)-2,5-二氢吡咯并[3,4-C]吡咯-1,4-二酮	2,5-Bis(2-ethylhexyl)-3,6-bis(5"-hexyl[2,2':5',2"-terthiophen]-5-yl)-2,5-dihydropyrrolo[3,4-c]pyrrole-1,4-dione	1093468-95-1
2,5-双(5-叔丁基-2-苯并恶唑基)噻吩	2,5-Bis(5-tert-butyl-2-benzoxazolyl)thiophene	7128-64-5
2,5-双(1-萘基)-1,3,4-恶二唑	2,5-Di(Naphthalen-1-Yl)-1,3,4-Oxadiazole	905-62-4
3,6-二噻吩基-吡咯并吡咯二酮	2,5-Dihydro-3,6-di-2-Thienyl-Pyrrolo[3,4-c]Pyrrole-1,4-Dione	850583-75-4
2,7-双[N-(萘基)苯氨基]-9,9'-螺二[9H-茛]	2,7-Bis[N-(1-Naphthyl)Anilino]-9,9'-Spirobi[9H-Fluorene]	932739-76-9
2,7-双[N,N-双(4-甲氧基苯基)氨基]-9,9'-螺二[9H-茛]	2,7-Bis[N,N-bis(4-methoxyphenyl)amino]-9,9-spirobi[9H-fluorene]	1138220-69-5
2,9-双(3-(二甲氨基)丙基)蒽[2,1,9-def:6,5,10-d'e'f']二异喹啉-1,3,8,10(2H,9H)-四酮	2,9-Bis(3-(Dimethylamino)Propyl)Anthra[2,1,9-Def:6,5,10-D'E'F']Diisoquinoline-1,3,8,10(2H,9H)-Tetraone	117901-97-0
2,9-双[3-(二甲氧化氨基)丙基]蒽[2,1,9-def:6,5,10-d'e'f']二异喹啉-1,3,8,10(2H,9H)-四酮	2,9-Bis[3-(Dimethyloxidoamino)Propyl]-Anthra[2,1,9-Def:6,5,10-d'e'f']Diisoquinoline-	1558023-86-1
2,9-二(萘-2-基)-4,7-二苯基-1,10-菲咯啉	2,9-Di(Naphthalen-2-Yl)-4,7-Diphenyl-1,10-Phenanthroline	1174006-43-9
N,N'-二甲基-3,4,9,10-茛四甲酰二亚胺	2,9-Dimethylantra[2,1,9-Def:6,5,10-D'E'F']Diisoquinoline-1,3,8,10(2H,9H)-Tetraone	5521-31-3
2-[5-[[4-[4-(2,2-二苯乙烯基)苯基]-1,2,3,3a,4,8b-六氢环戊[b]吡啶-7-基]亚甲基]-4-氧-2-硫氧代噻唑-3-基]乙酸	2-[5-[[4-[4-(2,2-Diphenylvinyl)Phenyl]-1,2,3,3a,4,8b-Hexahydrocyclopenta[b]Indol-7-yl]Methylene]-4-Oxo-2-Thioxothiazolidin-3-	652145-28-3
9-乙基咔唑-3-甲醛-二苯胺	3-(Diphenylhydrazonomethyl)-9-Ethylcarbazole	73276-70-7
3,3'-二[1,4]苯并恶唑并[2,3,4-kl]吩恶嗪	3,3'-Bi[1,4]Benzoxazino[2,3,4-Kl]Phenoxazine	1395881-58-9
3,3-b)二噻吩[3-氟-2[(2-乙基己基)羰基]噻吩并[3,4-b]噻吩二基)	3,3-B)Dithiophene}[3-Fluoro-2[(2-Ethylhexyl)Carbonyl]Thieno[3,4-	1380582-98-8
3,3-双[二(对甲苯基)氨基]联苯	3,3'-Bis[di(p-tolyl)amino]biphenyl	161485-60-5
3,4,9,10-四甲酰二亚胺	3,4,9,10-Perylenetetra-carboxylic Diimide	81-33-4
3,5-二(1-茛基)吡啶	3,5-Di(1-Pyrenyl)Pyridine	1246467-58-2
3,9-双(2-亚甲基-(3-(1,1-二氰基甲基)-6,7-二氟)-茛酮)-5,5,11,11-四(4-己基苯基)-二噻吩并[2,3-d:2,3-d]-s-茛并[1,2-b:5,6-b]二噻吩	3,9-bis(2-methylene-((3-(1,1-dicyanomethylene)-6,7-difluoro)-indanone))-5,5,11,11-tetrakis(4-hexylphenyl)-dithieno[2,3-d:2,3-d]-s-indaceno[1,2-b:5,6-b]dithiophene	2097998-59-7
4,4',4''-Tris[9,9-dimethylfluoren-2-yl(phenyl)amino]triphenylamine	4,4',4''-Tris[9,9-dimethylfluoren-2-yl(phenyl)amino]triphenylamine	303111-06-0
4,4'-双(4,6-二苯基-1,3,5-三嗪-2-基)-1,1'-联苯	4,4'-Bis(4,6-Diphenyl-1,3,5-Triazin-2-Yl)-1,1'-Biphenyl	266349-83-1
4,4-双(5-甲基-2-苯并噻基)二苯乙烯	4,4'-Bis(5-methyl-2-benzoxazolyl)stilbene	2397-00-4

中文名称	英文名称	CAS
4,6-双(3,5-二(吡啶-4-基)苯基)-2-甲基嘧啶	4,6-Bis(3,5-di(pyridin-4-yl)phenyl)-2-methylpyrimidine	1030380-51-8
4,7-双(5-正辛基-2-噻吩基)-2,1,3-苯并噻二唑	4,7-Bis(5-Octyl-2-Thienyl)-2,1,3-Benzothiadiazole	1171974-28-9
4,7-二(2-噻吩基)-2,1,3-苯并噻二唑	4,7-Di(2-Thienyl)-2,1,3-Benzothiadiazole	165190-76-1
EH-IDTBR	4-Thiazolidinone, 5,5'-[[4,4,9,9-Tetrakis(2-Ethylhexyl)-4,9-Dihydro-S-Indaceno[1,2-B:5,6-B']Dithiophene-2,7-Diyl]Bis(2,1,3-Benzothiadiazole-7,4-Diylmethylidyne)]Bis[3-	2102510-60-9
5''-([1,1'-联苯]-4-基)-1,1':4',1'':3'',1''':4''',1''''-对五联苯	5''-([1,1'-Biphenyl]-4-Yl)-1,1':4',1'':3'',1''':4''',1''''-Quinquephenyl	6326-64-3
5-氮杂螺[4.4]壬烷双(三氟甲磺酰)亚胺	5-Azoniaspiro[4.4]nonane Bis(trifluoromethanesulfonyl)imide	706759-44-6
5-Azoniaspiro[4.4]nonane Chloride	5-Azoniaspiro[4.4]nonane Chloride	98997-63-8
5-Azoniaspiro[4.4]nonane Iodide	5-Azoniaspiro[4.4]nonane Iodide	45650-35-9
5-氮杂螺[4.4]壬烷四氟硼酸盐	5-Azoniaspiro[4.4]Nonane Tetrafluoroborate	129211-47-8
7,7-二[1,4]苯并恶嗪并[2,3,4-kl]吩恶嗪	7,7'-Bi[1,4]benzoxazino[2,3,4-kl]phenoxazine	1395881-55-6
9,9-二[4-[二(萘-2-基)氨基]苯基]芴	9,9-Bis[4-(N,N-Bis-Naphthalen-2-yl-Amino)Phenyl]-9H-Fluorene	916061-87-5
(4-溴苄基)三苯基溴化磷	9,9-Bis[4-[N-(1-naphthyl)anilino]phenyl]fluorene	510775-24-3
9,9-二甲基-2,7-双[N-(间甲苯基)苯胺]芴	9,9-Dimethyl-2,7-bis[N-(m-tolyl)anilino]fluorene	143886-11-7
9-乙基咪唑-3-甲醛-N-苄基-N-苯胺	9-Ethyl-3-(N-Benzyl-N-Phenylhydrazonomethyl)Carbazole	75238-79-8
双(10-羟基苯并[h]喹啉)铍	Bis(10-Hydroxybenzo[h]Quinolinato)Beryllium	148896-39-3
碘化铋	Bismuth(III) Iodide	7787-64-6
C106 Dye	C106 Dye	1152310-69-4
溴化铯	Cesium Bromide	7787-69-1
氯化铯	Cesium Chloride	7647-17-8
碘化铯	Cesium Iodide	7789-17-5
苯并咪唑花	D-(+)-Melibiose Monohydrate	79534-91-1
二环己基(1-甲基-2,2-二苯基环丙基)膦	Dicyclohexyl(1-Methyl-2,2-Diphenylcyclopropyl)Phosphine	1240372-42-2
二茛并[1,2,3-cd:1',2',3'-lm]芘	Diindeno[1,2,3-Cd:1',2',3'-Lm]Perylene	188-94-3
二喹喔啉并[2,3-a:2',3'-c]吩嗪	Diquinoxalino[2,3-A:2',3'-C]Phenazine	214-83-5
十二烷-1-胺氢溴酸盐	Dodecan-1-Amine Hydrobromide	26204-55-7
硫氰酸胍	Guanidine Thiocyanate	593-84-0
己胺氢溴酸盐	Hexylamine Hydrobromide	7334-95-4
聚乙撑二氧噻吩-聚(苯乙烯磺酸盐)	high-conductivity grade	
试剂IEICO-2F	IEICO-2F	2055812-53-6
试剂ITIC-Th	ITIC-Th	1889344-13-1
试剂J51	J51	1393529-03-7
乙酸铅(II)三水合物	Lead Acetate Trihydrate(II)	6080-56-4
二氯化铅	Lead Dichloride	7758-95-4
碘化铅	Lead(II) Iodide	10101-63-0
2,2':6',2''-三联吡啶-4'-甲酸甲酯	Methyl 2,2':6',2''-Terpyridine-4'-Carboxylate	247058-06-6
甲胺三氟乙酸盐	Methylamine Trifluoroacetate	107326-29-4
甲基醋酸胺	Methylammonium Acetate	6998-30-7
N,N-双(2,5-二叔丁基苯基)-3,4,9,10-二甲酰亚胺	N,N-Bis(2,5-Di-T-BU-Phenyl)-3,4,9,10-perylenedicarboximide	83054-80-2
4,4-双[二(3,5-二甲苯基)氨基]-4-苯基三苯胺	N-(4-BiPhenyl)-N,N-Bis[4-[di(3,5-Xylyl)Amino]Phenyl]Amine	249609-49-2

中文名称	英文名称	CAS
N-(9,9-二甲基-9H-芴-2-基)-9,9-二甲基-N-苯基-9H-芴-2-胺	N-(9,9-Dimethyl-9H-Fluoren-2-Yl)-9,9-Dimethyl-N-Phenyl-9H-Fluoren-2-Amine	165320-27-4
N,N,N,N-四(4-甲基苯基)联苯胺	N,N,N',N'-Tetrakis(4-methylphenyl)benzidine	76185-65-4
N,N,N',N'-四([1,1'-联苯]-4-基)[1,1':4',1''-三联苯]-4,4''-二胺	N,N,N',N'-Tetra([1,1'-Biphenyl]-4-Yl)[1,1':4',1''-Terphenyl]-4,4''-Diamine	952431-34-4
N,N'-双(2,6-二异丙基苯基)-3,4,9,10-花四甲酰二亚胺	N,N'-Bis(2,6-Diisopropylphenyl)-3,4,9,10-Perylenetetracarboxylic Diimide	82953-57-9
N,N'-双(4-甲氧基苯基)-3,4,9,10-花四甲酰二亚胺	N,N'-Bis(4-Methoxyphenyl)-3,4,9,10-Perylenetetracarboxylic Diimide	6424-77-7
N,N'-二(4-联苯基)-N,N'-二苯基联苯胺	N,N'-Di(4-Biphenyl)-N,N'-Diphenylbenzidine	134008-76-7
N,N'-二正辛烷基-3,4,9,10-花四甲酰二亚胺	N,N'-Di-N-Octyl-3,4,9,10-Perylenetetracarboxylic Diimide	78151-58-3
N,N'-二戊基-3,4,9,10-花二甲酰亚胺	N,N'-Dipentyl-3,4,9,10-perylenedicarboximide	76372-75-3
N,N'-二(4-甲氧基-2-甲基苯基)-N,N'-二苯基联苯胺	N,N'-Bis(4-Methoxy-2-Methylphenyl)-N,N'-Diphenylbenzidine	169685-34-1
N,N-双(9,9-二甲基-9H-Fluoren-2-Yl)-N,N-二苯基联苯胺	N,N-Bis(9,9-Dimethyl-9H-Fluoren-2-Yl)-N,N-Diphenylbenzidine	361486-60-4
N,N-双[4-(二苯基氨基)苯基]-N,N-二苯基联苯胺	N,N'-Bis[4-(diphenylamino)phenyl]-N,N'-diphenylbenzidine	209980-53-0
N,N-二(1-萘基)-N,N,9,9-四苯基-9H-芴-2,7-二胺	N,N-Di(1-Naphthyl)-N,N,9,9-Tetraphenyl-9H-Fluorene-2,7-Diamine	357645-40-0
N,N'-二(2-萘基)-N,N'-二苯基-1,4-苯二胺	N,N'-Di(2-Naphthyl)-N,N'-Diphenyl-1,4-Phenylenediamine	139994-47-1
N1,N1'-([1,1'-联苯]-4,4'-二基)双(N1-(萘-1-基)-N4,N4-二苯基苯-1,4-二胺)	N1,N1'-([1,1'-Biphenyl]-4,4'-Diyl)Bis(N1-(Naphthalen-1-Yl)-N4,N4-Diphenylbenzene-1,4-Diamine)	910058-11-6
N1,N1-双(4-(二苯基氨基)苯基)-N4,N4-二苯基苯-1,4-二胺	N1,N1-Bis(4-(Diphenylamino)Phenyl)-N4,N4-Diphenylbenzene-1,4-Diamine	105389-36-4
N,N-二乙基-1,3-丙二胺二盐酸盐	N1,N1-Diethylpropane-1,3-Diamine	99310-71-1
N1,N4-二苯基-N1,N4-二间甲苯基苯-1,4-二胺	N1,N4-Diphenyl-N1,N4-Di-M-Tolylbenzene-1,4-Diamine	80223-29-6
N2,N2,N2',N2',N7,N7,N7',N7'-八对甲苯基-9,9'-螺联[芴]-2,2',7,7'-四胺	N2,N2,N2',N2',N7,N7,N7',N7'-Octa-P-Tolyl-9,9'-Spiro[Fluorene]-2,2',7,7'-Tetraamine	515834-67-0
N2,N7-二苯基-N2,N7-Di-M-Tolyl-9,9'-螺二[芴]-2,7-二胺	N2,N7-Diphenyl-N2,N7-Di-M-Tolyl-9,9'-Spiro[Fluorene]-2,7-Diamine	1033035-83-4
N,N,N',N'-四联苯基联苯二胺	N4,N4,N4',N4'-Tetra([1,1'-Biphenyl]-4-Yl)-[1,1'-Biphenyl]-4,4'-Diamine	164724-35-0
N4,N4,N4',N4'-四(4-甲氧基苯基)-[1,1'-联苯]-4,4''-二胺	N4,N4,N4',N4'-Tetrakis(4-Methoxyphenyl)-[1,1'-Biphenyl]-4,4'-Diamine	122738-21-0
N4,N4'-([1,1'-联苯]-4,4'-二酰基)双(N4,N4',N4'-三苯基-[1,1'-联苯]-4,4'-二胺)	N4,N4'-([1,1'-Biphenyl]-4,4'-Diyl)Bis(N4,N4',N4'-Triphenyl-[1,1'-Biphenyl]-4,4'-Diamine)	167218-46-4
N4,N4'-二(萘-1-基)-N4,N4'-二(萘-2-基)-[1,1'-联苯]-4,4'-二胺	N4,N4'-Di(Naphthalen-1-Yl)-N4,N4'-Di(Naphthalen-2-Yl)-[1,1'-Biphenyl]-4,4'-Diamine	374592-88-8
N4,N4-二([1,1'-联苯]-4-基)-N4'-萘-1-基)-N4'-苯基-[1,1'-联苯]-4,4'-二胺	N4,N4-Di([1,1'-Biphenyl]-4-Yl)-N4'-(Naphthalen-1-Yl)-N4'-Phenyl-[1,1'-Biphenyl]-4,4'-Diamine	897671-42-0
N,N'-二(4-甲基苯基)-N,N'-联苯二胺	N4,N4'-Diphenyl-N4,N4'-Di-P-Tolyl-[1,1'-Biphenyl]-4,4'-Diamine	20441-06-9
萘-1-基甲胺氢碘酸盐	Naphthalen-1-Ylmetanamine Hydroiodide	256222-11-4
正辛基碘化铵	N-Octylammonium Iodide	60734-63-6
PBDB-T	PBDB-T/PCE12	1415929-80-4
聚[2,5-双(3-十二烷基噻吩-2-基)噻吩并[3,2-b]	PBTTT-C12	888491-18-7



中文名称	英文名称	CAS
聚[[9-(1-辛基壬基)-9H-咪唑-2,7-二基]-2,5-噻吩二基-2,1,3-苯并噻二唑-4,7-二基-2,5-噻吩二基]	PCDTBT	958261-50-2
试剂PCE-10	PCE-10	1469791-66-9
聚(3,4-乙烯二氧噻吩)-聚(苯乙烯磺酸酯)	PEDOT:PSS	155090-83-8
并五苯	Pentacene	135-48-8
茈-3,4,9,10-四羧酸二酐	Perylene-3,4,9,10-Tetracarboxylic Dianhydride	128-69-8
聚[(5,6-二氟-2,1,3-苯并噻二唑-4,7-二基)-alt-(3,3''-二(2-辛基十二烷基)-2,2',5',2'',5'',2''''-四噻吩	PffBT4T-2OD	1644164-62-4
菲并[1,2-b:8,7-b]二噻吩	Phenanthro[1,2-b:8,7-b']dithiophene	1491133-64-2
苯乙溴化铵	Phenethylammonium Bromide	53916-94-2
苯硫氰酸铵	Phenylammonium Thiocyanate	27248-14-2
苯基双(4-(1-苯基-1H-苯并[d]咪唑-2-基)苯基)氧化膦	Phenylbis(4-(1-Phenyl-1H-Benzo[D]Imidazol-2-Yl)Phenyl)Phosphine Oxide	1426143-77-2
颜料红149	Pigment Red 149	4948-15-6
Platinum phthalocyanine	Platinum phthalocyanine	14075-08-2
聚[2,2'-[(2,5-双(2-辛基十二烷基)-3,6-二氧代-2,3,5,6-四氢吡咯并[3,4-c]吡咯-1,4-二基)]二噻吩-5,5'-二基-alt-2,2'-联噻吩-5,5'-二基]	Poly[2,2'-[(2,5-Bis(2-Octyldodecyl)-3,6-Dioxo-2,3,5,6-Tetrahydropyrrolo[3,4-C]Pyrrole-1,4-Diyl)]Dithiophene-5,5'-Diyl-Alt-2,2'-Bithiophene-	1267540-03-3
聚乙撑二氧噻吩	POLY(3,4-ETHYLENEDIOXYTHIOPHENE)	126213-51-2
四甲基丙烯酸酯封端的聚乙撑二氧噻吩	Poly(3,4-ethylenedioxythiophene), tetramethacrylate end-capped solution	
聚(乙撑二氧噻吩)-block-聚(乙二醇)溶液	Poly(3,4-ethylenedioxythiophene)-block-poly(ethylene glycol) solution	
聚(3,4-乙撑二氧噻吩)-聚(苯乙烯磺酸盐)	Poly(3,4-ethylenedioxythiophene)-poly(styrenesulfonate)	
3-己基取代聚噻吩	Poly(3-Hexylthiophene)	104934-50-1
聚(3-己基噻吩-2,5-二基)	Poly(3-hexylthiophene-2,5-diyl) (regioregular) [for organic electronics]	110134-47-9
聚(3-辛基噻吩-2,5-二基-co-3-癸氧基噻吩-2,5-二基)	Poly(3-octylthiophene-2,5-diyl-co-3-decyloxythiophene-2,5-diyl)	
聚(9,9-二辛基-9H-芴-2,7-二基)	Poly(9,9-Dioctyl-9H-Fluorene-2,7-Diyl)	195456-48-5
试剂Poly(benzimidazobenzophenanthroline)	Poly(benzimidazobenzophenanthroline)	34398-57-7
Poly(Fe-btpyb) Purple	Poly(Fe-btpyb) Purple	2288769-59-3
聚[[1,2,3,6,7,8-六氢-2,7-双(2-辛基十二烷基)-1,3,6,8-二氧代苯并[lmn][3,8]菲咯啉-4,9-二基][2,2'-联噻吩]-5,5'-二基]	Poly[[1,2,3,6,7,8-Hexahydro-2,7-Bis(2-Octyldodecyl)-1,3,6,8-Dioxobenzo[Lmn][3,8]Phenanthroline-4,9-	1100243-40-0
聚[(2,6-(4,8-双(5-(2-乙基己基-3-氟)噻吩-2-基)-苯并[1,2-B:4,5-B']二噻吩))-ALT-(5,5-(1,3-二-2-噻吩-5,7-双(2-乙基己基)苯并[1,2-C:4,5-C]二噻吩-4,8-二酮)	Poly[[4,8-Bis[5-(2-Ethylhexyl)-4-Fluoro-2-Thienyl]Benzo[1,2-B:4,5-B']Dithiophene-2,6-Diyl]-2,5-Thiophenediyl[5,7-Bis(2-Ethylhexyl)-4,8-Dioxo-4H,8H-Benzo[1,2-C:4,5-C]Dithiophene-1,3-Diyl]-2,5-Thiophenediyl]	1802013-83-7
聚[2,1,3-苯并噻二唑-4,7-二基[4,4-双(2-乙基己基)-4H-环戊并[2,1-B:3,4-B]二噻吩-2,6-二基]]	Poly[2,1,3-Benzothiadiazole-4,7-Diyl][4,4-Bis(2-Ethylhexyl)-4H-Cyclopenta[2,1-b:3,4-b]Dithiophene-2,6-Diyl]]	920515-34-0
聚[2,1,3-苯并噻二唑-4,7-二基-2,5-噻吩二基(9,9-二辛基-9H-芴-2,7-二基)-2,5-噻吩二基]	Poly[2,1,3-Benzothiadiazole-4,7-Diyl-2,5-Thiophenediyl(9,9-Dioctyl-9H-Fluorene-2,7-Diyl)-2,5-Thiophenediyl]	782469-77-6
聚[2,5-双(3-十四烷基噻吩-2-基)噻吩[3,2-b]噻吩	Poly[2,5-Bis(3-Tetradecylthiophen-2-Yl)Thieno[3,2-B]Thiophene	888491-19-8
聚苯胺	Polyaniline	25233-30-1
聚苯胺(翠绿亚胺盐)	Polyaniline (emeraldine salt)	

中文名称	英文名称	CAS
IEICO-4F	Propanedinitrile, 2,2'-[[4,4,9,9-Tetrakis(4-Hexylphenyl)-4,9-Dihydro-S-Indaceno[1,2-B:5,6-B']Dithiophene-2,7-Diyl]Bis[[4-[(2-Ethylhexyl)Oxy]-5,2-Thiophenediyl]Methylidene(5,6-Difluoro-3-Oxo-	2089044-02-8
氯化亚锡	Stannous Chloride	7772-99-8,(稳定提供多瓶与
Tetrabenzoporphine	Tetrabenzoporphine	52952-31-5
硫氰酸与十二烷基胺的化合物	Thiocyanic Acid Compound With Dodecylamine	22031-31-8
Tin II 2,3-Naphthalocyanine	Tin II 2,3-Naphthalocyanine	110479-58-8
碘化锡(II)	Tin(II) Iodide	10294-70-9
二氢基二氯化锡	Tin(IV)Phthalocyaninedichloride	18253-54-8
三(4-联苯基)胺	Tri(Biphenyl-4-yl)Amine	6543-20-0
三(五氟苯基)硼烷	Tris(Pentafluorophenyl)Borane	1109-15-5
三[4-(2-噻吩基)-4-联苯基]胺	Tris[4'-(2-thienyl)-4-biphenyl]amine	1092356-36-9
三[4-(2-噻吩基)苯基]胺	Tris[4-(2-Thienyl)Phenyl]Amine	142807-63-4
试剂Y6	Y6	2304444-49-1
β-丙氨酸氢碘酸盐	β-Alanine Hydroiodide	2096495-59-7

## 更多产品服务请咨询



业务一部

电话: 13817370081

QQ: 1952207695



业务二部

电话: 13817370080

QQ: 2402195779



业务三部

电话: 13817678381

QQ: 2574946592



业务四部

电话: 18301893121

QQ: 2102259500



业务五部

电话: 15821261143

QQ: 3119412031



业务六部

电话: 13310163795

QQ: 2316674096



业务七部

电话: 18301893120

QQ: 2689503877



市场部

电话: 13917439331

QQ: 1850926438